P7 HR&BP sports bracelet manual

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Manual English (original manual is with visual images)

1. Main menu:

Display with time / date / week

2. Heart rhythm menu:

Go to this menu by pressing the button when you feel a small vibration the green light goes on, the data appears after 15 seconds. The best way to measure your heart rhythm is to wear the bracelet 1-2cm behind the lump of your wrist. Make sure that the sensor is as close to the sheet as possible. If you use this function throughout the day, the data about your heart rhythm is kept on the application in the form of a graph.

Remember: a one-off heart rate test is not kept by the application or mobile phone.

3. Blood pressure menu:

Press the button until you feel a vibration with green light, the data will appear approximately 25 seconds later (wait patiently for the data to appear). If this function is used all day in combination with the application, the application keeps track of the data by means of a histogram.

Remember: One-time blood pressure test will not be kept in the application or mobile phone.

4. Pedometer: Also called pedometer and keeps track of how many steps are taken.

5. Distance: Keep track of how many km has been walked.

6. Calories: Keep track of how many calories have been burned.

7. Sportmode:

Press the button for some time to get to sports mode, if you are in sports mode press the button softly to change sport. For example, go running, swimming, basketball or badminton. When you are in a certain mode, press the button for some time to start the test. He starts after the small vibration. Under this mode various things are kept, such as heart rhythm, time, steps and calories burn. This depends on the sport you are doing. Press the arrow back to get out.

8. Turn the wrist for menu:

The screen gives light when the wrist is turned, if you want this function this can be enabled in the application.

9. Sleep test:

The bracelet automatically tests sleep, the application stores the sleep rhythm in this can be seen in the data history on the app.

10. Close menu

Go to the menu which is shown in the original manual, now you have to hold the button long enough until it closes.

Bluetooth connections

1. Download the application
2. You can download the application on your mobile (HPLUS)
3. The QR code can also be scanned on the manual or on the back of the box to download it
4. Open the Application choose the correct bluetooth name HPLUS and connect to the application. After a successful connection the smartband automatically synchronizes time, date and week of the mobile phone

Application specifications and functions

1. Pedometer (pedometer) (see visual images in the original manual):

To achieve your goals you can see the number of kilometers you have walked, the number of calories you have burned, the total number of steps you have taken and the application can also draw a conclusion about what needs to be improved.

2. Sleep (see visual images in the original manual):

The Smartband automatically tests the sleep, the sleep data can be read in the application. For example, the application analyzes how deep sleep is, this can vary in deep or light sleep. Then the application immediately gives a conclusion and advice.

3. Heart rhythm (see visual images in the original manual):

Press the button until you feel a vibration and see green light, it shows data after 15 seconds. You can set the heart rhythm that the application keeps a whole day, this is reflected in a graph in the application.

4. Blood pressure (see visual images in the original manual):

Press the button until you feel a vibration and see the blood pressure icon on the screen, data is displayed after 25 seconds. Here too you can configure the application to keep it for a whole day. The data of this can be seen again in a histogram in the application.

Set up:

My device: click on connect device to set the following functions.

- (a) Alarm: set a reminder for your alarm clock.
- (b) Reminder: set the time when you want to receive a reminder each time.
- (c) Screen saver time: Set when you want the screen to light up.
- (d) Message reminder: set a reminder function to incoming messages.
- (e) Synchronize data: click this function to synchronize data from the smartband.
- (f) Smartband unbundle: disconnect application and smartband.
- (g) Delete data from the smartband: delete all test data from the application.
- (h) Delete data from the application: delete the set data from the application.
- (i) Upgrade: click to download a new version.

Set goals:

- (a) Pedometer goal: set goals for each day.
- (b) Sleep goals: set sleep targets according to the personal situation.
- (c) Constitutional index: obtain the index by height and weight.
- (b) Fat burning heart rhythm: shows the effects of sport by analyzing heart rhythm.
- (c) Software version: shows the version of the application.