

Aphraheals

Sea Salt Band-Aid
to treat Infections



Heal your inflammation quickly
with proven miracle cure Aphraheals™

Contents of the box:



4 Bags with sea salt flour (Sea Salt & Fleur de Sel)

Our Celtic Sea Salt & Fleur de Sel comes from the Guérande region in Brittany, France and both forms are very rich in minerals and trace elements. The main minerals are magnesium, calcium, potassium and iron.



4 Band-Aids

To apply the gauze with saline solution to the skin.



1 pipette

To moisten the side of the Aphraheals that will cover the wound.



1 headband

To keep the Band-Aid and gauze with the Aphraheals solution in place while sleeping.

After years of research, we cracked the code.



Finally, no more allergic reactions and slow, complicated healing processes. There is a special ingredient in Aphraheals: the flower of sea salt. This is 15 times stronger than regular sea salt.

This unique method guarantees quick results.

How to make and use Aphraheals - Piercing Aftercare solution

Method 1:

Make an Aphraheals - solution and heal your sensitive piercing in 5 steps. For infected piercings see Method 2.

1. Wash your hands well with soap, remove any residue and dry your hands.
2. Open the Aphraheals - sachet
3. Place the sachet in a glass of clean, lukewarm water. After 10-15 minutes the content of the sachet will dissolve in the water. Dispose of the empty sachet.
4. Soak a cotton ball, piece of gauze or paper towel in the solution and place it on the affected area. Leave it on there for a minimum of 15 minutes. Leave it on for a longer period of time if your piercing still hurts after the recommended time.
5. Repeat this twice a day until the piercing is no longer sensitive/red and/or swollen.

Method 2:

Make an Aphraheals - solution and heal your infected piercing in 6 steps.

1. Wash your hands well with soap, remove any residue and dry your hands.
2. Open the Aphraheals - sachet
3. Fill a glass with lukewarm water and fill the pipette.
4. Moisten the side of the Aphraheals that will be placed on the wound.
5. Place the Aphraheals on the wound and cover it with a band Aid and headband to secure that it stays in place.
6. Let it work overnight and moisten the Aphraheals regularly by using the pipette. Repeat these steps until the infection is healed.

BEFORE



24



48



Order today via Aphraheals.com

Aphraheals - Celtic Sea Salt & Fleur De Sel

A common problem with healing piercings is the formation of a bump (due to moisture buildup) near the piercing. This can be small and harmless, but it can sometimes grow bigger and cause pain.

These types of bumps are often caused by constant pressure on the piercing (for example while sleeping) or a weakened immune system.

After years of research and several experiments, the Piercings Works Team has finally found a solution to this common problem.

It is a gauze pad filled with Celtic Sea Salt & Fleur De Sel. The gauze must be moistened before it is placed over the bump to ensure that the moisture is drawn from the bump.

This simple saline solution, that may seem too simple at first, is very effective! The gauze needs to be left on the infected area for at least a few hours.

The best result is achieved by letting it sit for an average of 8 hours, for example while sleeping.

In most cases, results can be seen immediately. Usually, the bump will have disappeared after 1 to 3 treatments with Aphraheals sea salt.

This solution works well with a titanium piercing rod or when you wear a NoPull piercing disc on the back of the piercing.

Adopt a healthy lifestyle, get sufficient sleep, and try to apply as little pressure on the piercing as possible.

About Aphraheals Celtic Salt & Fleur de Sel

Our Celtic Sea Salt & Fleur de Sel comes from the Guérande region in Brittany, France and both forms are very rich in minerals and trace elements. The main minerals are magnesium, calcium, potassium and iron. Of all salts, Celtic salt has the lowest sodium chloride content, making it the softest sea salt in the world. This can be attributed to the over two thousand year old salt extraction process.

The salt is extracted manually using a traditional method that is more than 2000 years old. Water from the cold sea current of the Atlantic Ocean is channeled through an ingenious channel system to shallow clay marshes.

The acidic clay in these swamps neutralizes the very basic seawater to an extent. Celtic sea salt with a sodium chloride content of about 90 percent is one of the most sodium-poor sea salt in the world.

The clay marshes also have a purifying effect. The seawater is channeled from one swamp to another, with the swamps acting as a natural purification method, because microbes in the clay neutralize impurities in the seawater. The entire process takes about two weeks, after which the salt water is collected in a basin where it evaporates by the sun, leaving only the salt crystals.

The salt is not washed, dried or refined, so the mineral balance remains completely intact. This purification method makes it 100 percent pollution-free.

Fast and simple healing for your piercing

It's great to have a beautiful, healed piercing, but what should you do when your piercing gets infected? Learn about 7 things you need to know to help your piercing heal smoothly.

1. A new piercing

The skin around the newly created piercing is often a little red and / or a little thicker than usual. It is also usually sensitive to the touch. That is part of the game and you do not have to and cannot do anything about it.

2. Simple and fast healing of your piercing

For a piercing to heal quickly and properly, personal and physical hygiene is of relevance. The area around the piercing needs to be thoroughly cleaned daily. It is important that all soap residue and shampoo are removed while showering as this can cause irritation and lead to infections.

3. Cleaning inaccessible places

To clean inaccessible places, it is useful to immerse your piercing in sea salt solution for 15 minutes. We recommend that you place a piece of gauze in the AphraHeals solution and place the moistened gauze on the piercing. You can't do this often enough. It is necessary to dry your piercing afterwards. If it stays wet much longer, it could lead to infections as well, as bacteria can build up. Swimming in a clean sea also helps!

4. Immune system

In addition, it is very important that you work on maintaining a strong immune system. Make sure to focus on the following:

- Quality sleep
- Limit stress
- Stick to a healthy and varied diet

5. Don't touch the piercing

It's very important to touch the piercing as little as possible. Do not put pressure on the piercing and don't sleep on it.

6. Professional piercing artist

All of the above will work as long as the piercing is the right size and shape. This is why it is important to have it done by a qualified piercer in a professional studio.

7. Proper materials

Finally, check if your piercing is made of the right material.

The appropriate materials for a piercing are titanium, surgical steel 316L, solid gold and all other hypo-allergic and highly polished materials.

If you follow these recommendations during the healing process of your piercing, you minimize the chance of infection. If, after applying all these tips, your piercing is still sensitive, I would recommend a special solution: AphraHeals.

We have been using this miracle cure for 10 years and it works. We don't call it "the infection killer" for no reason! Usually, an infected piercing will heal after a day when you using AphraHeals. The swelling will go down, there will no longer be any discharge, and the color of the skin around the piercing will go back to normal. Curious to find out more and experience the power of AphraHeals yourself?

10 ways to use Aphraheals

Below you'll find the 10 uses of AphraHeals and how it can help you treat various issues and stay healthy in times of COVID.

1. Treat infected piercings

AphraHeals is a perfect solution for infected piercings and provides a quick fix now that piercing studios are closed and GPs are difficult to reach.

2. High magnesium content

AphraHeals is rich in magnesium. Magnesium strengthens your immune system and reduces stress.

3. Lasts up to a week

If you dissolve AphraHeals in a bottle of water, you can use it for a period of 7 days to care for your piercing.

4. You can eat it

Did you run out of salt? Aphraheals consists of the highest quality of sea salt flour, and can therefore also be eaten.

5. Top restaurants use it

Top restaurants often use the flower of sea salt instead of regular sea salt for their meals because of its unique taste.

6. Household disinfectant

Use the AphraHeals solution as a disinfectant for your home.

7. Wound healing

AphraHeals helps to heal all types of wounds.

8. Helps with blemishes

AphraHeals also helps with blemishes.

9. Can be used as a scrub

It can be used as a scrub by mixing it with coconut or baby oil.

10. Mouthwash and toothpaste

It can be used as a mouthwash by dissolving the sachet in a glass of water. Additionally, it can be used as a toothpaste. To make your own toothpaste, mix AphraHeals with baking soda / sodium bicarbonate.



BEFORE



AFTER



Order today via Aphraheals.com

Why do people no longer use piercing sprays after using Aphraheals?

When we started piercing 20 years ago, things were very different. There were far fewer piercings and the most popular piercing was the navel piercing. We pierced a lot of people and mostly young people came to our studios or found us at big festivals. At the time, piercings were not so widely accepted in society. It was rebellious, daring, unusual. We loved it.

Severe inflammation or infection

In most cases, piercings healed beautifully and with proper care and attention, the healing process went smoothly and quickly. Unfortunately this wasn't always the case. We had acquired a lot of experience and were curious about what went wrong. Sometimes we would encounter bad infections or overgrown piercings. How could this be?

We always used clean materials and had a lot of experience. We were the first studio to be officially approved by the GGD and followed all hygienic guidelines. The GGD advised to use antibacterial soap for healing. So we gave this to our customers. This wasn't a success... For some, it was hard to properly rinse off the soap which led to infections. It's safe to say that this wasn't an ideal solution.

We continued to think of alternatives. A customer told us that she applied Tea Tree oil to her piercing. She noticed a positive effect and her piercing seemed to heal quickly. Great idea! Or was it? Yes, for most people this could provide some relief but others experienced severe allergic reactions. We still hadn't found a solution that would work for all.

Dissolving sea salt

We continued to try out different things and found that dissolving sea salt in water helped and caused no allergic reactions. This seemed to be going in the right direction! Proper healing unfortunately still took quite some time so there was definitely room for improvement.

We went on to investigate sea salt further. We found out that there were different types of sea salt and made different types of solutions using Celtic sea salt. This worked better, but still not fast enough. Then we discovered the so-called “flower of sea salt”. It contains 15 times more minerals than ordinary sea salt. After a long period of experimenting, we discovered an effective way to use this sea salt.

We gave this sea salt to people with severe inflammation and overgrown piercings and advised them to use it in the same way. We could not believe what happened! This miracle cure caused fast healing. Sometimes results were achieved within just 1 night of use! Sometimes it took a little longer, but the healing process was faster than we had ever seen.

This was a real breakthrough!

We were really happy with the results. Customers of our Amsterdam studio who experienced infections told others about this remedy and we were able to help more and more people effectively. And because we want to help even more people spend up their healing process, we have decided to make this product available for sale. This is how AphraHeals was born! The miracle cure for fast and effective piercing healing.

How can you use AphraHeals?

It's simple: If you have an infection, moisten the bag filled with the "flower of sea salt" and place it on the infected area. Keep it in place with a band-aid or hair band (or both). Leave it on overnight.

When the wound dries, add a little water using the included pipette. After a few nights, your infection will be a thing of the past. If you don't have inflammation and you want to speed up the healing process you can simply use AphraHeals for 15 minutes and repeat this twice a day. This way the piercing will heal quickly and smoothly.



Frequently Asked Questions about Aphraheals

How long does an Aphraheals (AH) bag last?

You can use the bag in two different ways depending on which method you choose.

Method one: Are you using AH to treat and infected piercing? Then we recommend that you use it once and throw it out afterwards.

Method two: Are you using AH as aftercare for a piercing that is not infected? You can make a solution using AH and water and use this solution for a period of 7 days.

Can AH cause an allergic reaction?

This rarely happens. The bag itself is made of Abaca, a banana variety from the Philippines. The plant is also known as Manila hemp. In addition, the most important ingredient of AH is the flower of sea salt and hardly anyone is allergic to this.

How did Aphraheals come about?

We took years to develop AH in our Amsterdam studio, Piercings Works.

Can I also treat pimples using AH?

You certainly can! This remedy is particularly good for acne. You will achieve the best results when you leave the Band-Aid with AH on all night. You can use it in the same way as you would with an infected piercing. Usually the pimple dries out and becomes pus-free after one night.

Is it normal to experience itchiness in the areas where I use AH?

Yes, it is normal to experience different sensations in the area where you apply an AH patch. These sensations indicate activity, extracting pus from your wound, and calming any inflammation.

MEDICAL SILICONE PIERCING DISC



Are you bothered with a piercing that gets too far in your skin while you have a rod that is absolutely long enough?

Medical Silicone Piercing Discs offer a solution!

These discs prevent your piercing from getting too far in your skin, but that's not all! Medical Silicone Piercing Discs also work perfectly with cartilage bumps. Silicone helps reduce scar tissue while reducing the irritation of the jewelry.

Are you a fan of Aphraheals?

Are you fan of Aphraheals: aftercare for your piercing?

Apply to become an ambassador, promote Aphraheals and get paid!

Sounds good, right? Start immediately:

1. Sign up; it's free
2. You will instantly receive further instructions by email
3. Receive an unique link, only for you
4. Send your unique link with your friends on Facebook or in your own blog. Promote as often as possible!
5. Wait ... and watch sales come in!

Sign up here and start promoting:

aphraheals.com/become-an-ambassador/

KALAFEERO



1.

With clean hands, submerge the pipette into a clean glass of (lukewarm) water and squeeze to fill up with water.

.....



2.

Remove the bag from the seal pack and squeeze the pipette to transfer 3ml of water.

.....



3.

Cover the affected area with the wet part of the bag. Place the band-aid on top and peel off the sides to make it stick to your skin.

.....



4.

Carefully squeeze water under the band-aid using the pipette when it's drying up. For affected areas around the head, use the headband if necessary.

Repeat this process until the wound is healed