

PAX Rescue and Transport Systems Vacuum





it's all about performance







Your dependable partner in the field

We only use specially selected and manufactured materials for our PAX Vacuum Products, always according to our specifications.

In combination with our product philosophy "it's all about performance", we developed products with a high the field.

Discover your PAX Vacuum Product:

PAX Vacuum Products General Information	Page 4-5
PAX Vacuum Splints Notes and product variants	Page 6-11
PAX Vacuum Mattresses Notes and product variants	Page 12 - 21
PAX Vacuum Spare Parts Helpful Accessories	Page 22 - 29
PAX Vacuum Products Special Special products	Page 30 - 35



PAX Products for daily use

You must be able to rely on your PAX Product, just as the patient must be able to rely on your skills.

For over 20 years, we at PAX have been working on the quality of our products with our high standards of design and functionality using innovative and resilient materials. The use of new technologies as well as the expertise gained from practical experience and the close dialogue with our customers are used to develop a PAX Product. The design, colour selection and the use of materials adapted to the intended purpose result in a PAX Product that is designed to meet the highest demands. We're only satistied if the PAX Product stands the test and you as customer are pleased with it as well.

In this product brochure "PAX Rescue and Transport Systems" we present a variety of our new and revised products. You can find more innovative PAX Products at www.pax-bags.com.





degree of utilisation for rescuers as well as the patient. Our lightweight and robust materials in combination with our design and the cross-system components make our PAX Vacuum Products your dependable partner in



Layered construction

Our PAX Vacuum Mattresses are characterised by a multi-chamber system with a 2-layer structure. The EPS beads of the layer on the patient side can be spread over a large area and thus enable the mattress to be shaped effectively to the patient. The chamber system of the layer remote from the patient is devided in a way that guarantees the necessary stability of the moulded vacuum mattress.



Image simila







As compact as you need it

In addition, the layered design also ensures that you can stow everything away in the smallest spaces.

Well packed is half the battle

Our transport bags are especially designed to fit the PAX Vacuum Products and thus facilitate space-saving storage and transport. With them, your PAX Vacuum Product is always well protected and hygienically packed. They are available in 3 different sizes.

and even radiolucent.

All PAX Vacuum Products are tested according to EN 1865-1 and have very good X-ray permeability.

PAX Vacuum Products

Practical - well thought out

Feels good

This significantly increases patient comfort.



Because it has to be clean

Our materials are tested for intensive cleaning and disinfection. You can rely on that. The carrying handles are no exception. They have closed pores and are therefore perfect for cleaning. This also applies to the belt holders which are welded in from the inside to maintain a smooth surface.

Our PAX Vacuum Products should be checked for damage and cleaned again after each use. After cleaning them with a damp cloth, disinfect the PAX Vacuum Products with a standard disinfectant.

PAX Vacuum Products are practical



The use of our PAX-Flex Tec MP material ensures that the surface of your PAX Vacuum Mattress or PAX Vacuum Splints is soft and supple.



Simple and safe - with push-in fasteners

Our PAX Vacuum Products do not use a Velcro system for belt fixation system and use the PAX Buckles instead - this is easier to clean. Especially when opening, this prevents the unintentional movement of the patient's extremities. The well-being of the patient is in the centre of our attention.



PAX Vacuum Splint -

Forearm



PAX Vacuum Splint

- Rec - Arm

PAX Vacuum Splint -Elbow



PAX Vacuum Splint -Leg



PAX Vacuum Splint -Rec - Joint



PAX Vacuum Splint -Rec - Leg







PAX Vacuum Splints

Our system for immobilisation of the extremities

Our two different PAX Vacuum Splint sys- Coming in two different material combinatems enable a fast and effective immo- tions and shapes, the splints are very useful bilisation of the extremities in daily use. aids for many situations.







PAX Vacuum Splint - Forearm

Our easy-to-fit PAX Vaccuum Splint - Forearm is fastened securely to the patient with exchangeable straps and is made of the soft and comfortable PAX-Flex Tec MP material.

51x34 **■** 0,35



Item no. 155545210

PAX Vacuum Splint - Leg

With our easy-to-fit leg splint, the patient can be treated safely. By closing the lower buckles, the PAX Vacuum Splint - Leg is preformed 3-dimensionally and thus achieves very good stabilisation. The splint is securely fastened to the patient with the interchangeable fixation belts. The PAX Vacuum Splint - Leg is made of the very soft and comfortable PAX-Flex Tec MP material.



plate for the individual labelling

C€ MD

FT

This product series is characterised

by its particularly low weight!

PRO +

回旅潟回

PAX Vacuum Splint - Elbow

by connecting the lower clasps before application. This enables effective application and shaping of the splint. The splint is securely fastened to the patient by means of the exchangeable fixation straps. The PAX Vacuum Splint - Elbow is made of our soft and comfortable PAX-Flex Tec MP material.



Item no. 157500000S

PAX Vacuum Splints - Set

Of course, our PAX Vacuum Splints are also available as a set. This allows you to store all designs in our space-saving PAX Transport Bag which is already included in the delivery.

PAX Vacuum Splints - Set consisting of: 1 x PAX Vacuum Splint - Leg 1x PAX Vacuum Splint - Forearm 1x PAX Vacuum Splint - Elbow 1 x PAX Transport Bag





PAX Vacuum Splint - Rec - Arm

For quick application to the arm, the PAX Vacuum Splint - Rec -Arm has a simple shape and is fastened to the patient by means of fixation straps with buckles. The underside of the splint consists of our strong PAX-Tec Plus MP, while the patient side is made of the soft and adaptable PAX-Flex Tec MP material. This facilitates easy adaptation to the patient's anatomy.

47x57 0.60



C€ MD

FT

PRO +



Item no. 162035210

PAX Vacuum Splint - Rec - Joint

The PAX Vacuum Splint - Rec - Joint consists of the robust PAX-Tec Plus MP material on the underside and the soft and adaptable PAX-Flex Tec MP material on the patient side. The simple shape in combination with the fixation straps including buckles facilitates quick application to the patient.





Be flexible in the field - whatever size you're using

You never know what to expect when you're in action, that's why our PAX Vacuum Mattresses are designed to accommodate a wide range of patient sizes and weight classes. This is ensured, among other things, by the material of the layered and multi-chamber system.

A class of its own - a good feeling

Our PAX Claim you can rely on.

-

Sophisticated arrangement for better usability

Your PAX Vacuum Mattress has an ingenious arrangement of the buckles. They are on the sides and not in the middle of the mattress. This means that the immobilised patient does not suffer pressure due to buckles lying on the mattress.

5 is trumps

PAX Fixation Straps have a colour code corresponding to the buckles. This helps you to "keep a cool head" and minimises the risk of misuse in terms of integrated safety for patients and users.

One for all

In the PAX Vacuum System of our PAX Vacuum Mattresses, only one belt system is used. You do not need 3 or 4 different systems, but can always

use one belt (also from the accessories) for all components in the system.

Included in scope of delivery: This belt is the universal belt extension that can be used individually.



PAX Vacuum Mattresses

Quality in form

The advantages and product develop- patients with suspected bone and spinal tion. Our PAX Vacuum Mattresses fulfil the ments of our PAX Vacuum Splints are column injuries as well as cranialbrain or recommendations of the S3 guidelines continued in our PAX Vacuum Mattresses. polytrauma. They envelop the patient and polytrauma. They are suitable for the immobilisation of provide very good support and stabilisa-





for the individual labelling



Item no. 155165210

PAX Vacuum Mattress - I-Mat - Handles

The PAX Vacuum Mattress - I-Mat - Handles is our entry-level product and is manufactured with the same PAX Care as all our products. The straight cut is reminiscent of an uppercase "I" and optimised for use in the medical sector. The PAX Vacuum Mattress - I-Mat - Handles is supplied with 8 loop-in handles and our 5 colourcoded fixation straps. The strap and carrying loop eyelets are recessed for better cleaning and the carrying handles are completely closed-pored. To avoid pressure points on the patient, the buckles of the PAX Fixation straps are placed on the outer sides of the mattress. The individually usable strap extension is also included in the scope of delivery of this model. Our PAX Vacuum Mattress - I-Mat - Handles fulfils the requirements of the EN 1865-1 standard as well as the recommendations of the S3 guideline - polytrauma and is made of the very soft and comfortable PAX-Flex Tec MP material.

70x202 4,23

 $\stackrel{1x}{=} 20x65x95 \stackrel{2x}{=} 27x65x60$







Ω

Ω

Π

0

 \cap

Customised labelling card

You decide what you want to label

PAX Vacuum Mattresses offer a clear field for your individual identification and labelling - well visible, large enough and transparently protected. We offer you a format template for individual labelling, which you can design on our website.



Product-specific . labelling card





ltem no. 155145210

PAX Vacuum Mattress - Mummy-Mat - Handle Bar

With the PAX Vacuum Mattress - Mummy-Mat - Handle Bar, we orientated ourselves on the shape of the body and shaped it as an oversized "V". The "V" oversize enables a better securing of the patient due to the larger span from the pelvis. This version has a handle bar for transport. This vacuum mattress can be moulded around the patient beyond the sides. The patient is secured with the help of the colour-coded fixation straps. The strap and carrying loop eyelets are recessed for even better cleaning and the carrying handles are completely closed-pored. To avoid pressure points on the patient, the buckles of the PAX Fixation straps are placed on the outer sides of the EN 1865-1 standard as well as the recommendations of the S3 guideline - polytrauma. The side facing the patient is made of the very soft and comfortable PAX-Flex Tec MP material. The matching transport bag and the individually usable belt extension are already included in the scope of delivery.

97x202 4,70



ltem no. 155155210

PAX Vacuum Mattress - Mummy-Mat - Handles

This version of the PAX Vacuum Mattress - Mummy-Mat - Handles is identical in construction with the version `handle bar' but has loop-in handles for transport. With the PAX Vacuum Mattress - Mummy-Mat - Handles we orientated ourselves on the body shape and shaped it as an oversized "V". The "V" oversize enables a better securing of the patient due to the larger span from the pelvis. It can be modelled around the patient beyond the sides, the patient is secured with the help of the colour-coded fixation straps or fixed in place. For cleaning and hygienic reprocessing, the strap and carrying loop eyelets are are completely closed-pored. To avoid pressure points on the patient, the buckles of the PAX Fixation straps are placed on the outer sides of the mattress. An individually usable belt extension as well as the matching transport bag are already included in the scope of delivery.

97x202 4,70

1x 15x90x95 23x90x62







Just in case: PAX Repair Patch Vacuum Item .no. 157510000



PAX Vacuum Mattress - Ergo-Mat - Handle Bar

The PAX Vacuum Mattress - Ergo-Mat - Handle Bar is the flagship of our new PAX Vacuum Mattresses. In addition to the anatomical V-cut, it has an attached head section which allows the use of the vehicle stretchers patient straps without any problems. It has colourcoded fixation straps with high viz yellow buckles, which are colour-coded to the straps on the PAX Logo, as well as an individually usable belt extension. The strap and carrying loop eyelets are completely closed-pored for better cleaning. To avoid pressure points on the patient, the buckles of the PAX Fixation straps are placed on the outer sides of the mattress. Our PAX Vacuum Mattress - Ergo-Mat - Handle Bar fulfils the requirements of the EN 1865-1 standard and the recommendations of the S3 guideline - Polytrauma and is made of the very soft and comfortable PAX-Flex Tec MP material on the side facing the patient.

97x203 4,80

$\stackrel{1x}{=} 20x90x97 \stackrel{2x}{=} 23x90x60$



Item no. 155135210

PAX Vacuum Mattress - Ergo-Mat - Handles

The PAX Vacuum Mattress - Ergo-Mat - Handles is identical in construction with the PAX Vacuum Mattress - Ergo-Mat - handle bar, but has, as the name suggests, loop-in handles instead of the handle bar for transport. It complies with the recommendations of the S3 guidelines and is made of the very soft and comfortable PAX-Flex Tec MP material at the side facing the patient. The colour-coded fixation straps are used for quick identification. For optimal cleaning and hygienic reprocessing, the strap and carrying loop eyelets are completely closed-pored. The PAX Vacuum Mattress - Ergo-Mat - Handles has, in addition to the anatomical V-cut, an attached head section which enables the use of the vehicle stretchers belts without any problems.









///



PAX Vacuum Mattress - Ergo-Mat - Handle Bar & Head Support

For this model "PAX Ergo-Mat with integrated head fixation and neck support", we attached a handle bar for transport. Our new, patented head-fixation system and the integrated neck support prevent the hyper-reclination of the cervical spine during transport. The head fixation system is vacuumed together with the mattress and thus gives the head considerably more support. It also has colour-coded fixation straps for quick colour-coded fixation straps with high-viz yellow buckles. The strap and carrying loops are recessed for even better cleaning and the carrying handles are completely closed-pored. An individually usable belt extension, the matching transport bag and the fixation strips are already included in the scope of delivery. The fixing strips are also available as a set of 10 spare parts.

97x203 5,20



Item no. 160875210

PAX Vacuum Mattress - Ergo-Mat - Handles & Head Fixation

For this model "PAX Ergo-Mat with integrated head fixation and neck support" we attached handles for transport. Our new, patented head fixation system and the integrated neck support prevent hyper-reclination of the cervical spine during transport. The head fixation is suctioned together with the mattress and thus gives the head considerably more support. It also has colour-coded fixation straps with high-viz yellow buckles for quick identification. The belt and carrying loop eyelets are recessed for even better cleaning and the carrying handles are completely closed-pored. The fixing strips are also available as a available as a set of 10 spare parts.

97x203 5.20







Compatible for all extraction systems

The PAX Valve is supplied in two versions. Version for use of the valve with a fingertip: with this a piece of tubing is attached to the valve, the inner diameter of which allows the use of standard fingertips. Variant for use of the valve with the PAX Vacuum Pump: This variant allows the hose of the PAX Vacuum Pump to be plugged in directly without an additional adapter.

PAX Vacuum Spare Parts Useful accessories

tem belt fasteners, belt extensions, special hardly anything to be desired.

For the PAX Vacuum Products you get a hose connections, carrying handles, repair lot of accessories, which are coordinated stickers, transport bags and the highly effiwith the entire system. With its cross-sys- cient PAX Vacuum Pump the system leaves

Air out - without any effort

PAX Vacuum Pumps are made of high-quality and robust materials. Through the good sliding properties of the piston in the aluminium cylinder and the short stroke distance the PAX Vacuum Pumps enable a fast and effective suction of your PAX Vacuum Products. The PAX Vacuum Pump is available in two versions.

PAX Vacuum Pump - Hand: For the use of the small and easy to store vacuum pump.

PAX Vacuum Pump - Foot: This pump is equipped with an additional foot bar made of robust stainless steel. This enables effective suction with only one hand.

> CE MD

When the air is out...

...it should stay outside.

If a mishap should occur during use and the PAX Vacuum Products draw air, it is not an unsolvable problem. Our self-adhesive, transparent and easy-to-use PAX Repair Patch Vacuum can be affixed to the leak in a matter of seconds and your PAX Vacuum Product will immediately work again as usual.









Item no. 157510000

PAX Repair Patch Vacuum

Self-adhesive PAX Repair Patch for the simple repair of minor damage to a PAX Vacuum Product. For repairs on the patient side as well as on the underside of our PAX Vacuum Products.

10x7.5



Item no. 160881508

PAX - Head Fixation PAX Vacuum Mattress

The PAX Replacement Velcro Strips of the PAX Head Fixation for our PAX Vacuum Mattresses consisting of: 1x 25er Velcro strip high-viz yellow, length 50 cm, 1x 50er Fleece strip high-viz yellow, length 50 cm, with oval cut-out in the middle, (2 x 12,5 cm). Available in accessories as a set of 10.

0.02



Item no. 164244203

PAX Adapter T-Piece for Head Fixation

PAX Adapter T-piece for mounting the head fixation on the PAX Vacuum Mattress. Made of 50 mm webbing with hook strap as T-piece. 1 set consists of 2 individual T-pieces each.

0,35





03

PAX Vacuum Valve

Views with PAX Safety Loop



PAX Vacuum Valve

Views with PAX Safety Loop directly on the product







Item no. 274270000

PAX Vacuum Valve with Adapter for Fingertip

Consisting of valve for PAX vacuum products with glued-in PVC tube for attaching a fingertip. Supplied without safety loop and fingertip. 0,03





PAX Vacuum Valve Mattress and Splint

Spare part PAX Vacuum Valve for PAX Vacuum Mattress and PAX Vacuum Splint.

0.44x0.26 0.07



Item no. 268790303

PAX Safety Loop

Vacuum Valve

PAX Safety Loop for fastening of the PAX Vacuum Screw Valve Mattress and Splint to the PAX Vacuum Products. Available in a set of 10.



5,2x11,4 0,01





for hand suction.

10x30,7x5,9 0,70

it's all about performance

Item no. 273130003



CEMD PAX Vacuum Pump - Foot 03 Even with the foot, our PAX Vacuum Products are quickly ready for use. With the help of the PAX Vacuum Pump - Foot, they can be attached to the patient in the twinkling of an eye. 1,13 C€ MD PAX Vacuum Pump - Hand 03 With the help of the PAX Vacuum Pump - Hand, our PAX Vacuum Products can be modelled around the patient. Here the version Item no. 156990008 PAX Vacuum Mattress Handle Set (10 pcs.) Set consisting of 10 single handles. Replacement handles for easy assembly to PAX Vacuum Mattresses. Proven and ergono-

mic PAX Handle.

Item no. 162420008

PAX Vacuum Mattress Handle Set (2 pcs.)

CEMD

Replacement handles for easy mounting on PAX Vacuum Mattresses. Proven and ergonomic PAX Handle. Set consisting of 2 handles.

Item no. 157490410

PAX Transport Bag Vacuum Splints

The entire PAX Vacuum Splint Set can be stored in this matching transport and storage bag and thus be conveniently transported.

35x50x14 0.34

Item no. 157470410

PAX Transport Bag Vacuum Mattress

Of course we also offer you the matching transport and storage bags for our PAX Vacuum Mattresses. Here they find sufficient space and can be transported comfortably.

71x77x14 1,60



Find more information on our website as well as new products in the shop and more information in the download area. Let yourself be inspired.

Do you still have questions?

You can reach our consultants by phone or you can arrange an appointment with one of our field employees.

We are there for you - Your reliable partner in everyday life and in the field.









The PAX Fixation Straps of the vacuum mattresses and the leg and arm splints are fastened to the patient by means of straps and buckles to hold them in place and in shape. and in shape. The removable straps can easily be cleaned and disinfected or exchanged.

Item no. 157004200

PAX Fixing Strap Set Vacuum Mattresses

Fixation straps for easy and quick mounting on the PAX Vacuum Mattresses and easy fixation of the patient. Unique colour code for quick identification even in difficult situation or poor lighting conditions. Colour-coded buckle in high-viz yellow.

4,5 1,40





Included in scope of delivery: This strap is the universal strap extension that can be used individually.

ltem no.	Article name	Colour	Dimension (cm)
277554210	PAX Fixing Strap Extension Grey - Vacuum Mattress		111
277544203	PAX Fixing Strap Black - Vacuum Mattress		147
277534232	PAX Fixing Strap Blue - Vacuum Mattress		167
277524204	PAX Fixing Strap Yellow - Vacuum Mattress		187
277514201	PAX Fixing Strap Red - Vacuum Mattress		187
277504211	PAX Fixing Strap Green - Vacuum Mattress		177

Item no. 274731110

PAX Fixation Strap Set - Vacuum

Splint Elbow

This replacement strap set consists of:

- 2 x 50 cm fixation belt with plug-in buckle and PAX-Flex Tec 3 mm T-fastening piece
- 2 x 75 cm fixation belt with plug-in buckle and PAX-Flex Tec 3 mm H-fastening piece

tem	no.	274741110	

PAX Fixation Strap Set - Vacuum

Splint Forearm

This replacement strap set consists of: 3 x 75 cm fixation belt with plug-in buckle and PAX-Flex Tec 3 mm H-fastening piece



- and 2 ladder buckles for attachment
- and 2 ladder buckles for attachment





C E MD



Item no. 274751110

PAX Fixation Strap Set - Vacuum

Splint Leg

This replacement strap set consists of:



3 x 75 cm fixation belt with plug-in buckle and PAX-Flex Tec 3 mm T-fastening piece 1 x 50 cm fixation belt with plug-in buckle and PAX-Flex Tec 3 mm T-fastening piece





PAX Vacuum Products Special

Vacuum mattresses for special applications

We developed our PAX Vacuum Products well. We rise to the challenge and adapt as a standard range. But as you are used the product characteristics and specific reto from PAX, there are special designs as quirements to your purposes.





ican QR Code and watch the film bout our PAX Rescue Bag RTS All





PAX Rescue Bag - RTS AIR

The PAX Rescue Bag - RTS AIR is distinguished by the consistent omission of Velcro in order to protect the safety-relevant ropes and straps for winching. With the all-round zip, the PAX Rescue Bag - RTS AIR can be operated by one person, even in adverse conditions. The internal lacing system in the head area and the colour-coded compression straps allow the PAX Rescue Bag to be quickly adjusted to the patient.

The redundantly designed loop system for attaching the hanger to the helicopter for transport as an external load provides extra safety. Vertical rescue or abseiling is also possible and permitted. All loops and ropes are safely stowed in pockets. The PAX Rescue Bag - RTS AIR is tested according to the requirements of EUROCAE ED-14G / RTCA DO-160G for Loose Equipment and fulfils the "Category U2 - Unknown Helicopter Frequencies". The tests have been carried out according to DIN EN 60068-2-64:2009 (vibration test) and DIN EN 60068-2-27:2010 (shock test) in the test laboratory TÜV SÜD.

The PAX Rescue Bag - RTS AIR can therefore be used in flight operations without further verification.

information the product:

Scan QR Code

to get more

Item no. 155585210

PAX Vacuum Mattress - AR 2

We combined particularly light and hard-wearing materials to meet the requirements of air rescue. The multichamber 2-layer filling gives this mattress sufficient stiffness. Because of the frequently encountered ski boots, we also processed the outer material in the foot area of the patient side. The handles are fixed directly to the mattress to save weight and space. In the neck area, an insert supports the particularly sensitive area.



The PAX Vacuum Mattress - AR 2 is designed for the PAX Rescue Bag RTS AIR.







PAX Vacuum Mattress AR 1

The special challenge in the development of this vacuum mattress lay in the labelling and the maximum weight of the product. The disinfectant-resistant labelling in the form of a special foil, which we produced using a specially manufactured tool and a highfrequency welding, was only one of the challenges. Thanks to special materials we managed to reduce the total weight of the mattress to 3,5 kg despite the large weight. The customer uses this mattress in the rescue bag, on the helicopter stretcher and for transport. The fixation straps are made of single-use Velcro straps to save weight.

90x205 3,50

 $\stackrel{1x}{=}$ 16x90x97 $\stackrel{2x}{=}$ 20x90x60



Special design of the marking

Jsing a special tool, single-colour ogos and /or lettering can be eproduced on the vacuum mattres

With the high-frequency process the ogo is connected inseparably and esistant to disinfectants on the mat-ress with a special foil. This process is characterised by its pecial durability.

egal Ar

Item no. 277025210

PAX Vacuum Mattress - Mummy-Mat - Plus - Handles & Head Fixation

Within the scope of a nationwide standardisation, we developed this product according to customers specifications and approved it as a medical product. We met the special requirements in terms of weight and stowage volume by means of the PAX Materials used, the special compartment division and the patented head fixation.





10		
FT		
PRO		

CEMD

ltem no. 268905310

PAX Vacuum Support - Roll-up Stretcher LAS

The PAX Vacuum Support - Roll-up Stretcher LAS as a positioning aid for the London Ambulance Stretcher has an adapted shape, welded, indented handles and a shoulder belt cut-out. It has a low filling capacity, as it only serves as a storage aid.

69x195 **2**,10







ltem no. 270134510

PAX Transport Bag for LAS Vacuum Support 39x59x21 **■** 0,30



C€ MD

10





Scan QR Code and tretcher LAS:



PAX it's all about performance

In this brochure you will find a selection of our range. The item number given is supplemented by your choice of colour. Of course, you can also make this selection in our online shop at www.pax-bags.com.

General symbols

H×B×F

KG

1x____

2x

PRO SEG

⁺<u></u>

†**__**__

Colours:



Made in Germany

high-viz yellow

Product is compliant with the applicable EU directives

Product is a medical device and complies with EU Directive 2017/745 - MDR

Materials for lasting functionality



Film material stretchable in all directions with a pleasantly soft feel, which makes it easy to model at low and high temperatures. Disinfectable and biocompatible, non-toxic to cells according to DIN EN ISO 10993-5.



This material quality is used for light products and indoors. Resistant, waterproof and disinfectable. Weight per unit area between 150 and 200 g/m².



Very robust PES fabric with film lamination on both sides. The outer side has a stronger coating to increase abrasion resistance. Good disinfectability and biocompatibility, non-toxic to cells according to DIN EN ISO 10993-5.

can QR code and view more inform n the PAX Materials



PAX Cleaning Service

You want us to wash and disinfect your PAX Product?

For organisations and companies we offer the PAX Cleaning Service. You can also send us your backpacks and bags for cleaning. For the duration of the cleaning we offer appropriate rental products. You can find out more about this on our homepage:www.pax-bags.com/en/about-pax/pax-service/





PAX Quality

Development

We have our own team of developers who work successfully with the countless developers in our customer base. We are all united by the desire to develop a PAX Product in the best possible way for its intended use.

Manufacture

In our manufactory we realise individual pieces - even the most For more than 15 years, we have been working with our individual customer wishes - and test our ideas for feasibility. partner. PAX Quality is taken very seriously there and a lot of attention to detail is put into every PAX Product.

Warranty

If you have registered your PAX Products for the extended warranty, you can use the advantage of the registration for to identify your product and get the product data from our system.

Simply enter the serial number that is on your PAX Product. If you have registered for a free warranty extension (see the special regulations at www.pax-bags.com), you will receive an attractive loyalty bonus of 10% on all services that relate to registered products within the scope of the repair services. For the duration of the repair, we will gladly provide you with a replacement product. The request for this replacement product must be stated within the scope of the repair registration. We will then get in touch with you.

Small batch production

In our production in Poland, our colleagues have been working with us for many years and bring the proven PAX quality into the PAX Products.

Large-scale production





The serial number required for registration can be found on the handle of many of our products or clearly visible inside the product.

PAX Repair Service

Sometimes even our products suffer the scuffs of everyday use. And sometimes something gets broken. No problem for you, because the PAX Repair Service is at your side.





And this is how the **PAX Repair Service works:**

Simply register in our PAX Online Shop under "My PAX Repair Service" in your personal account area to use our repair tool. Here you can choose whether you want us to repair your product or have it professionally cleaned and disinfected.

In the next step, you enter the PAX Product to be repaired / cleaned by entering its item no or the product name in the search field and confirm your selection on the following page. If you have decided on a repair you can now select the repair reason(s) before the price is displayed and additional services such as loan products and pick-up service are offered. The same applies if you have decided on a cleaning service, except that here the repair reasons are omitted.

On the last page, the shopping basket, the services you have selected are clearly summarised. By clicking on "Proceed to checkout", you complete your order bindingly.

If this area is not yet visible to you, please contact us, we will then activate it for you.









These product brochures are valid from 09/2021. All product photos and colours shown in them are for visualisation purposes only and are subject to colour variations due to printing technology. They are not binding in terms of colour. All products are delivered without contents. Some of the fillings shown are for information purposes only. The image material comes from various sources (shutterstock.com). The brochures have been prepared to the best of our knowledge. Any errors or technical changes due to technical





reasons, which serve progress, cannot be ruled out. presented by the managing directors: Thomas Busch, Local Court Oldenburg HRA 3860; the GmbH is re- 30345339.7).

www.pax-bags.com

https://de-de.facebook.com/paxbags/

https://www.instagram.com/paxbags/

The technical values are laboratory values and have Andreas Harms - Sales tax identification number acno claim to correctness. The current warranty condi- cording to § 27a turnover tax law: DE 220516085 tions as well as the currently valid general terms and - Tax number: FA Oldenburg 64/206/02304 - Reconditions can be found on our homepage www. gister court: Local court of Oldenburg - Commercial pax-bags.com Published by: X-CEN-TEK GmbH & register number: HRA 3860 - PAX is a trademark Co. KG - Westerburger Weg 30 - 26203 Warden- registered at the German Patent and registered burg X-CEN-TEK GmbH & Co. KG is represented trademark at the German Patent and Trademark by the personally liable partner: X-CENTEK GmbH, Office in Munich (registration number/ file number:



X-CEN-TEK GmbH & Co. KG • Westerburger Weg 30 • D-26203 Wardenburg Phone: +49 (0)44 07 - 7 14 76 - 0 • Email: info@pax-bags.de • www.pax-bags.com/en

it's all about performance