

# The mask that protects **YOU.**

## Protection Mask



On the job, at home or out and about, we all have to think about virus protection. It's particularly important during a pandemic and wherever viral load is high. Masks are an essential component in hygiene rules. How well the mask protects depends on its fit and filtration.

The advantages of our Protection Mask:

- no-leak fit **protects you and others**
- **top filter performance**
  - filters out the tiniest particles ( $\varnothing 0.1 \mu\text{m}$ )
  - similar to an FFP2 mask
- **one size** fits (almost) all
- **reusable**  
(can be disinfected and sterilized)

The respected Nelson Labs tested the filter material.

- ✓ **Passed the virus test**
- ✓ **Passed the bacteria test**

**VIRUS  
PROTECTION  
CONFIRMED**

Benefit from Löwenstein Medical's decades of expertise in the development of respiratory masks.



### Reprocessing

Mask part	Material	Chemical Disinfection*	Thermal Disinfection**	Sterilization ***	Home environment Sterilize ****	Cycle
Mask cushion	Silicone	✓	✓	✓	✓	30
Headgear	Silicone	✓	✓	✓	✓	30
Fastening clip	Plastic Polyoxymethylene (POM)	✓	✓	✓	✓	30
Filter insert	Plastic Polyamide (PA)	✓	✓	✓	✓	30
Filter holder	Plastic Polyamide (PA)	✓	✓	✓	✓	30
Particle filter	Fleece	-	-	-	-	At the latest after one shift or when breathing resistance increases

\* Sekusept Aktiv 2%: 15 minutes, Gigasept FF (new) 7%: 15 minutes  
 \*\* 94 °C: 5 minutes  
 \*\*\* 134 °C: 7 minutes  
 \*\*\*\* 100 °C: 15 minutes

### Technical Data

Product class per EU Medical Devices Directive 2017/745	I
Dimensions (W x H x D)	89 mm x 105 mm x 70 mm
Weight	93 g
Temperature range: - Operation - Transport and Storage	+5 °C to + 40 °C -20 °C to +70 °C
Allowable humidity: Transport and Storage	5% to 95% relative humidity
Period of use	up to 8 hours (one shift) or change particle filter earlier if breathing resistance is excessively high
Service life - Mask - Particle filter (un-used)	30 cleaning cycles or 5 years after manufacture 5 years

#### Product description

Protection Mask  
Set of 10 particle filters

#### Article Nr.

LMT 26685  
LMT 26694



Löwenstein Medical Technology  
 Kronsaalsweg 40, 22525 Hamburg  
 Germany  
 T: +49 40 54702-0  
 F: +49 40 54702-461  
 info@loewensteinmedical.de  
 www.loewensteinmedical.de

Löwenstein Medical  
 Arzbacher Straße 80, 56130 Bad Ems  
 Germany  
 T: +49 2603 9600-0  
 F: +49 2603 9600-50  
 info@hul.de  
 www.hul.de

