**Posters scientifiques
housses anti-acariens**
***Scientific posters
anti-dust mite covers***

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**Impact of the use of dust mite covers manufactured with ISO 13485:2016 certified processus on the quality of life of allergic patients.*

A follow-up study carried out for 52 patients. C. Marti n, A. Vanlaeys, L. Salles, T. Porée - WAC 2019

Our Laboratory

OptimHal-ProtecSom is a research company that designs and manufactures in France medical devices, designed for patients suffering from respiratory diseases. The company, whose production unit, research center, and headquarters are based in Normandy, has become a global reference in the field of dust mite allergy protection and aerosol therapy. The company's success is the result of a team effort that brings together healthcare professionals, researchers, salespeople, administrative staff, technicians, and assembly workers.

Research as a key word

OptimHal-ProtecSom teams are constantly seeking innovations to create the devices of tomorrow. These devices will help to improve the quality of life for patients suffering from respiratory diseases. In this way, OptimHal-ProtecSom has developed numerous patented technologies in the field of inhalation.

Our various publications

In this catalog, you will find various publications on anti-dust mite covers. ProtecSom dust mite covers are a real textile barrier, impervious to allergens, uncoated and effective over time. They isolate patients from their dust mite allergens and reduce their symptoms by 62%*.

**Impact of the use of dust mite covers manufactured with ISO 13485:2016 certified processus on the quality of life of allergic patients.*

A follow-up study carried out for 52 patients. C. Marti n, A. Vanlaeys, L. Salles, T. Porée - WAC 2019

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Bibliographie

- Linneberg, A. et al. Burden of allergic respiratory disease: a systematic review. Clin Mol Allergy 14, (2016).
- Pingitore, G. & Pinter, E. Environmental interventions for mite-induced asthma: a journey between systematic reviews, contrasting evidence and clinical practice. Eur Ann Allergy Clin Immunol 45, 74-77 (2013).
- Juniper, E.F., Thompson, A.K., Ferrie, P.J., and Roberts, J.N. (2000). Development and validation of the mini Rhinconjunctivitis Quality of Life Questionnaire. Clin. Exp. Allergy J. Br. Soc. Allergy Clin. Immunol. 30, 132-140

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Poster CFA 2017

Evolution of rhinitis symptoms in mite-allergic patients using allergen-impermeable bedding covers: a multicenter, randomized, placebo-controlled, double-blind study

E. Furon ¹, F. Valingot-Anfray ², N. Boukhettala ¹, D. Carvaillo ³, A. Corbet ¹, N. Goldstein ⁴, D. Herbin ⁵, Y. Massabie ⁶, I. Mercier ⁷, T. Porée ¹, E. Puillandre ⁸, P. Rufin ⁷.

¹Laboratoire ProtecSom - Valognes (France), ²Allergologue - Bois Guillaume (France), ³Pneumologue Hôpital Privé de Chantilly - Chantilly (France), ⁴Allergologue - Avignon (France), ⁵Allergologue - Cherbourg (France), ⁶Allergologue - Marseille (France), ⁷Allergologue - Paris (France), ⁸Allergologue - La Teste-de-Buch (France).

Introduction

Allergic rhinitis is frequently under-diagnosed, despite its significant socio-economic impact¹. Patients generally remain symptomatic despite treatment. The first step in treating allergy is to reduce exposure to allergens. In dust mite allergy, although study results are inconsistent and clinical benefits controversial, the use of "anti-dust mite" mattress covers reduces exposure to allergens². The aim of this study is to provide objective data on the appropriateness of using integral woven cotton covers in patients with severe dust mite allergy and associated rhinitis.

Study method

Thirty mite-allergic patients were included in a multicenter, randomized, double-blind, placebo-controlled study. An allergen-impermeable mattress cover was installed on the beds of patients in the active group, and a non-impermeable cover on the beds of patients in the placebo group. The study consisted of: an inclusion consultation with clinical assessment, installation of the covers and a second clinical assessment after 12 weeks of uninterrupted use of the covers (figure 1). The severity of rhinitis was assessed using the Rhinoconjunctivitis Quality of Life Questionnaire (miniRQLQ)3 and the potential impact of pharmacological treatments was evaluated using the Medication Score Questionnaire (MSQ).

Bibliography

- Linneberg, A. et al. Burden of allergic respiratory disease: a systematic review. Clin Mol Allergy 14, (2016).
- Pingitore, G. & Pinter, E. Environmental interventions for mite-induced asthma: a journey between systematic reviews, contrasting evidence and clinical practice. Eur Ann Allergy Clin Immunol 45, 74-77 (2013).
- Juniper, E.F., Thompson, A.K., Ferris, P.J., and Roberts, J.N. (2000). Development and validation of the mini Rhinoconjunctivitis Quality of Life Questionnaire. Clin. Exp. Allergy J. Br. Soc. Allergy Clin. Immunol. 30, 132-140

Results

The miniRQLQ was completed by 25 patients (aged 12 to 68), and data from 22 patients could be analyzed (12 in the placebo group and 10 in the active group, figure 2). After 12 weeks, 9 out of 10 patients perceived a significant improvement in their symptoms (miniRQLQ delta score > 0.7) in the active group, versus 5 out of 12 in the placebo group (figure 3). A significant difference (table 1) in symptom improvement was observed between the active group (+1.53 ± 0.88) and the placebo group (+0.23 ± 1.19), particularly for nasal manifestations (2.3 ± 1.44 versus 0.25 ± 1.46). Medication scores remained similar in the 2 groups throughout the study (figure 4).

Conclusion

This study shows an improvement in the symptoms and quality of life of patients suffering from allergic rhinitis after using an anti-dust mite mattress cover made from woven cotton, and thus suggests a clinical benefit to their use.


Figure 1: Diagram of the study process.


Figure 2: Patient flow diagram.

	Placebo Group	Active Group
Δ total	0.23 ± 1.19	1.53 ± 0.88 *
Δ activities	0.11 ± 1.04	1.63 ± 0.96 **
Δ practical issues	0.29 ± 2.19	2 ± 1.87
Δ nose issues	0.25 ± 1.46	2.3 ± 1.44 **
Δ Eye issues	-0.39 ± 1.67	1.7 ± 1.82 *
Δ Other issues	0.92 ± 1.42	0.17 ± 0.88
Nb of patients Δ > 0.7	5	9
Nb of patients (Δ < -0.7)	2	0
Nb of stable patients	5	1


Figure 3: Averages of symptom scores per patient obtained with miniRQLQ0 et miniRQLQ13.


Figure 4: Drug score. Treatments, taken by participants throughout the study, were prioritized according to their ability to reduce allergic rhinitis symptoms: oral, conjunctival and nasal antihistamines = 1; nasal and inhaled corticosteroids = 2; oral corticosteroids = 3.

Control of the bedroom environment during the study


Figure 5: Complementary allergen avoidance measures monitored during the study using a follow-up questionnaire.

Environmental reasons for non-inclusion:

- Access to animals in the room
- Presence of plants in the room
- Use of box springs and/or headboards

Environmental reasons for exclusion:

- Change of bed for more than one week during the study period

Allergen avoidance measures implemented during the study (figure 5):

- Relative humidity
- Ventilation
- Regular washing of bedding
- Presence of « dust nests »

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Impact of processes certification on the liability of anti-dust mites bed covers

Thierry Porée, Nabile Boukhattala, PhD, Emeline Furon, PhD
R&D, Laboratoire Protec'som, Valognes, France

In relation to this presentation, I declare that there are no conflicts of interest.

Anti-dust mites covers

Simply killing dust mites does not remove allergens, several measures are essential to lower allergens and improve symptoms of allergic patients. Anti-mites covers are the most relevant device for allergen evicton¹.

Anti-mite barrier covers are manufactured from textiles with filtration properties usually validated by different certification labels. Although only textile certifications are mandatory, some covers are manufactured using processes ISO 13485² certified, subjecting covers to more stringent quality control measures.

The objective of the study was to evaluate the benefit of ISO 13485 certification on the liability of the final product.

Conformity assessment procedures

Manufacturers of class 1 medical devices are only required to provide a self-declaration of safety and performances, while ISO 13485 certification requires well-defined processes and strict recording (Figure 2).


Figure 2 : ISO 13485 full Quality assurance system

Methods

In this study, three different batches of micro-woven (MWC) and non-woven polyester polyamide (NWP) textiles went through the quality control processes required by the norm ISO 13485.

Textiles pore sizes were measured using an optical microscope (Bresser, USA).

Textile permeability was tested using a Rotomitest, apparatus composed of two compartments separated by the sample. Der p1 allergens were placed in compartment 1 and the apparatus was set to rotate for 18 hours (Figure 3).

Der p1 allergens that passed through the sample to compartment 2 were then measured with a Der p1 ELISA kit (Citeq Biologics, Netherlands).

Medical devices

In Europe, medical devices are designed and manufactured in compliance with directives 93/42/CEE³. Devices complying with those directives are entitled to display CE marking.

There are four classes of medical devices, ranging from low risk to high risk (Figure 1). Anti-dust mites covers are class one medical devices.


Figure 1 : European classification and regulatory requirements for medical devices.

ISO 13485 certification

International Organization for Standardization (ISO) standard 13485 presents the requirements for a quality management system for the design and manufacture of medical devices.

The norm ISO 13485 guarantees, with regular audit :

- The mandatory justification of the medical benefit.
- Product conformity to quality standard.
- Risk analysis throughout products development.
- Strict controls on all processes, from design to manufacturing.
- The control of the communication by a certification approved by European health authorities.

R&D@laboprotectsom.com

Results

Microphotographs of textiles (Figure 4 A,B and C) clearly show that anti-dust mites certified textiles (micro-woven/Texaal® Cotton and non-woven polyester polyamide/Noxaalon®) have a higher density of fiber when compared to cotton used to make common bed covers. Permeability testing of those fabrics (Figure 4 D) shows that both anti-dust mites certified textiles possess the required filtration properties.

Figure 3 : Rotomitest principle. Allergens are placed in compartment 1 and collected in compartment 2 after 18h.


Figure 4 : Representative microphotographs and filtration properties of tissues used. A. common cotton, B. Micro-woven cotton (MWC), C. non-woven polyester polyamide (NWP) and D. Permeability results of common cotton, MWC and NWP obtained using the rotomitest apparatus, showing the amount of allergen that passed through the fabrics, from compartment 1 to compartment 2.

Three batches of both anti-dust mites certified textiles were tested according to ISO 13485 monitoring processes. NWP fabrics showed irregular pore sizes, with 2 batches having pore size greater than 5 µm (Figure 5) while MWC showed similar pore sizes for all 3 batches. Permeability results (Figure 6) show that 2 batches of NWP fabrics are significantly more permeable to Der p1 allergen compared to MWC tissue (23,6ng +/- 0.7 vs 3,4ng +/- 0.02).


Figure 5 : Microphotographs of 3 batches of NWP fabrics and 3 batches of MWC fabrics, showing 2 batches of NWP with irregular pore sizes (> 5 µm).


Figure 6 : Permeability results showed that NWP fabrics are significantly more permeable to Der p1 allergen compared to MWC tissue (23,6ng +/- 0.7 vs 3,4ng +/- 0.02).

Thus, ISO 13485 quality controls allowed for dismissal of 2 out of 6 batches of unsatisfactory NWP textiles.

References

- ¹Lau S. Allergen Avoidance as Primary Prevention: Con." Clin Rev Allergy Immunol 2005 ; 28: 17-23.
- ²ISO 13485:2012 : Medical devices – Quality management systems..
- ³Council directive 93/42/EEC from the council of the European Communities of 14 June 1993 relating to the placing on the market of medical devices.

Conclusion

In addition to the initial test on textiles, as for medical devices of higher class, monitoring manufacturing processes certification is necessary to ensure the quality of the finished product, especially when using non-woven fabric. Therefore, ISO 13485 certification is a relevant criteria for anti-mite covers quality and thus effectiveness.

ISO 13485 certification is not required for anti-mite covers, but it is a pledge of quality

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Traduction du poster en français : "Impact of the use of dust mites covers manufactured with ISO 13485 certified processes on the quality of life of allergic patients - A follow-up study carried out on 52 patients"

Impact of the use of dust mites covers manufactured with ISO 13485 certified process on the quality life of allergic patients – A follow-up study carried out of 52 patients

C.Martin¹, A.Vanlaeys¹, L.Salles¹, T.Porée¹

¹Laboratoire OptimHal-ProtecSom - Valognes (France)

Background

- **Pathology:** allergic rhinitis¹ (prevalence: 31% of the French population)
- **Major cause:** dust mites
- **Symptoms:** breathing difficulties, nasal discharge and cough
- **Solution:** dust mites covers are tightly woven to create a mechanical barrier²
- **Purpose:** to evaluate the evolution of rhinitis symptoms of patients who used dust mites covers designed and manufactured according to the standards of the medical device (ISO 13485:2016)

Study material

Several brands of dust mites covers have been used by our patients (Texal Cotton, Texal Polyester, Noxaalon, and Noxaalon Bamboo). These covers are not chemically treated, and their anti-mite properties rely solely on mechanical filtration (mesh size: 0.5-5 µm).

Study method

- Fifty-two of the 107 questionnaires distributed to allergic people who have given their consent to attend the survey were returned.
- The questions were designed to assess patient satisfaction with the use of dust mite covers as well as the symptoms they experienced. The level of discomfort felt was assessed using a score of 0 (no discomfort) to 10 (severe discomfort).
- Values were expressed as mean ± standard error of the mean (SEM). Parametric Student test or non parametric Wilcoxon test was used to evaluate statistically differences depending on the distribution and the variance of samples (** p-value <0.001, * p-value <0.01, * p-value <0.05).

Results


Figure 1: Decrease of discomfort score from the first wave (between August and December 2018)


Figure 2: Decrease of discomfort score from the second wave (between January and August 2019)

The first wave (for 6 months) had determined that the level of discomfort felt by patients decreased considerably after the use of dust mites covers: 7.35 ± 0.2 vs 3.9 ± 0.2, i.e. a decrease in symptoms of 47%. The second wave of questionnaire (for 8 months) had determined that the level of discomfort felt by the patients decreased considerably after the use of dust mites covers: 7.58 ± 0.2 vs 2.88 ± 0.2, i.e. a decrease in symptoms of 62%. The difference in decreasing discomfort score between these two waves (47% VS 62%) is particularly explained by the season (autumn and winter, respectively)².


Figure 3: Evaluation of discomfort score per symptoms before and after implementation of dust mites covers

The discomfort score was systematically lower (from -58% for itchy skin to -72% for eye inflammation) in patients who used dust mites covers independently of the etiology of rhinitis.


Figure 4: Evaluation of discomfort score related to the moment of implementation of other elimination ways added to the use of dust mites covers

The implementation of other elimination ways seems to decrease the discomfort score (tendency). Moreover, no patient has implemented these evicton measures after the introduction of dust mites covers.


Figure 5: Distribution of allergic people who attend the survey according to age

Main users of dust mites covers are children less than 10 years old (26,32%) and people between 25 and 45 years old (27,63%).


Figure 6: Evaluation of discomfort score per age group before and after implementation of dust mites covers

The discomfort score was significantly lower after the implementation of dust mites cover independently of the patient age. However, the effect was more pronounced for children less than 10 years (-71%) and 25-45 years (-69%) groups.

Discussion

As all patients implemented other evicton measures, it is difficult to isolate the influence of dust mites covers from the influence of other evicton measures. Another study³ performed in "double-blind against placebo covers" has demonstrated the effectiveness of dust mites covers independently of other evicton measures. The follow-up study provides ongoing evaluation of the benefit of these measures and interesting information on the application of a set of preventive measures.

Conclusion

Using evicton measures, such as dust mites covers, enables the decrease of allergic rhinitis symptoms.

It would also be relevant to carry out a study distinguishing between patients those who have already implemented evicton measures and from those who have not yet done so, in order to quantify the decrease related solely to the use of ProtecSom dust mites covers.

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OptimHal-ProtecSom SAS