

Miracle Collar Success Stories!

"How would you like to inherit 57 horses overnight with most of them being cribbers? That's exactly what happened when we rescued a group of Thoroughbreds from a Davison, Michigan farm. If we hadn't gotten the Miracle Collar, the barn would have been gone and the fences right after!"

Barbara Baker
President of Horses' Haven

Horses' Haven is a non-profit organization located in Michigan that's dedicated to the welfare of horses and other large animals.



"...I have prescribed the Miracle Collar for all obsessive cribbers...the pressure to stop air-sucking can be applied without undue restriction of the throat area, thus removing a possible source of injury to this sensitive area. Let me thank you for providing a useful tool to treat obsessive cribbers."

Wayne Hachten, D.V.M.
East Penn Veterinary Associates, Inc.

"This wonderful invention is truly a miracle. I am so pleased with how well the collar quickly worked..."

Amanda Wenham
Wenatchee, WA

"...It is great that you have come up with such a remarkable product to help our wonderful equine companions. Thanks."

Prudence K. Heany
Two Time Tack
Bethesda, MD

"Thank you for creating what truly is a Miracle Collar. My Appaloosa, Spot, used to crib all the time, but since I have used your collar he doesn't even try. Spot wears it almost all the time — in his stall, turned out and to and from shows. I ride him English and sometimes I even leave the collar on him when riding. I recommended the Miracle Collar to my friends who have a very intense cribber named Chomps (who got his name from being such an intense cribber). All the other collars haven't worked on him, but I bet your collar will!"

Melissa Kurtz
Rochester, NY

We'd love to hear your success stories, questions or comments about the Miracle Collar, please call or write.

Toll Free in the U.S. and Canada: **800-932-8371**

Local and International: **330-674-1782**



7540 CR 201, PO Box 68
Mt. Hope, OH 44660-0068
www.ridethebrand.com

The One and Only Original Design

Miracle Collar®



Important Information About Adjusting and Caring for Your New Miracle Collar




100% SATISFACTION GUARANTEED

STOPS CRIBBING IMMEDIATELY!

Stop Cribbing Now With These Easy Steps:

The Miracle Collar must be adjusted correctly to stop your horse from cribbing. Please remember to check the fit frequently to ensure the collar's effectiveness.

- 1 Standing on the near side of your horse, place the wide bottom piece under the jaw with buckles facing you. The **CENTER TAB** should point forward between the jaws.
- 2 Position the shorter **NECK STRAP** over the neck, through the mane (or over the bridle path) and buckle when it is just snug enough to keep collar in place. You should always be able to slip a flat hand's width between the neck strap and your horse's neck. **NOTE: The neck strap does not apply the pressure to stop cribbing and should be adjusted only snug enough to keep the collar in place.**
- 3 Run the longer **BROWBAND STRAP** across the forehead just in front of the ears and buckle moderately tight. As when putting anything new on your horse, watch how your horse reacts to the Miracle Collar's adjustment and make further adjustments one hole at a time.



- 4 With your horse's head in the upright position, adjust the browband strap until it applies just enough pressure to stop your horse from cribbing.
- 5 To make sure the browband strap is adjusted tight enough not to drop over your horse's face, watch your horse eat off the ground. If the browband strap appears loose, continue to tighten it until it can't slide down over your horse's face and eyes. If necessary, punch additional holes in the strap for more adjustment.

Attention: Check the collar on a daily basis as leather becomes softer and more supple from daily use. The browband strap may need to be adjusted occasionally in the head down position described above. The browband strap should always fit tighter than the neck strap.

- 6 Check the neck strap again and make sure you can slip a flat hand's width between the neck strap and your horse's neck.
 - 7 Run straps through metal keepers on the buckles.
- Attention:** The browband strap applies the pressure to stop cribbing. Be sure it is properly adjusted at all times. The neck strap simply helps hold the Miracle Collar in place.

100% SATISFACTION GUARANTEED

We are deeply committed to providing you with high quality, dependable products that offer you years of use and satisfaction at a great value. And if any Weaver Leather product is not exactly what you expected or hoped for, simply return it promptly to the place of purchase for complete satisfaction. At Weaver Leather, we're not satisfied unless you are.

Caring for Your Miracle Collar

The Miracle Collar is constructed of premium quality harness leather. With proper care and maintenance, this collar will be effective and provide years of dependable service.

As a part of your daily grooming schedule, check the Miracle Collar for a proper fit as the leather may stretch and lengthen over time and need to be adjusted. Also inspect the Miracle Collar for any needed cleaning or repair. When the collar is fitted properly and kept clean, your horse should not develop sores. Clean and condition regularly to keep the Miracle Collar in shape. **A daily inspection of the Miracle Collar is absolutely necessary to help avoid any problems associated with the accumulation of dirt, manure, grime, sweat or problems that may occur from incorrect adjustment or position of the collar.**

Because horses are powerful and sometimes unpredictable animals, we recommend the Miracle Collar be used by experienced horse owners, following all directions and using common sense.

Miracle Collar[®]
WEAVER
LEATHER

Please keep these instructions for future references.



100% SATISFACTION GUARANTEED

ORIGINAL DESIGN

STOPS CRIBBING IMMEDIATELY!