CEKA Multi Set # 6700

2 x 4 g CEKA Multi Composite
6 automix tips
5 ml CEKA Multi Primer
Brush

CEKA Multi O Set # 6710

2 x 4 g CEKA Multi O Composite
6 automix tips
5 ml CEKA Multi Primer
Brush

CEKA Multi Composite # 6701

2 x 4 g CEKA Multi Composite
6 automix tips

CEKA Multi O Composite # 6711

2 x 4 g CEKA Multi O Composite
6 automix tips

CEKA Multi Primer # 6702

5 ml

CEKA Multi Mix Tips # 6705

20 pcs

BONDING STRENGTHS

Zirconium dioxide:	25.5 N/mm ²
Titanium:	21.4 N/mm ²
Gold:	21.7 N/mm ²
Cr-Co:	37.7 N/mm ²

AXIAL

ATTACHMENTS

CEKA AXIAL
PRECI-CLIX AXIAL
PRECI-CLIX RADICULAR
PRECI-BALL

EXTRACORONAL

ATTACHMENTS

CEKA EXTRACORONAL
PRECI-VERTIX
PRECI-CLIX EXTRACORONAL
PRECI-SAGIX
PRECI-S2

PROFILE

ATTACHMENTS

PRECI-BAR
PRECI-HORIX
PRECI-CLIP

DIGITAL

LIBRARIES

SPECIALTIES

PRECI-PROFILE
PRECI-POST
CEKA SOL
CEKA SITE
CEKA BOND
3C-BOND

CEKA Multi - CEKA Multi O

PERMA-RET
PRECI-SEP
PLASTICWAX
EXPANDO
MEASURING GAUGE

SIDE EFFECTS, WARNING AND PRECAUTIONS

- The attachments are intended for single use.
 - The products are non-sterile.
 - There is a risk of poor fit when patient conditions change.
 - Bacterial adhesion can be avoided by applying hygiene measures.
 - Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
 - The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
 - For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
 - Do not use products containing titanium.
 - Do not use products containing nickel in case of nickel allergy.
- The accessories RE H 79 and H 35 must be used outside the mouth.

CE0197

Prospectus no 297.M1-ed.E
© Copyright ALPHADENT NV 02/2019
BTW BE 0414.789.321 RPR Koninkrijk

CEKA PRECI-LINE

WWW.CKPL.EU

CEKA PRECI-LINE

SPECIALTIES

CEKA Multi / CEKA Multi O

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



ALPHADENT NV
Mannebeekstraat 33,
8790 Waregem, Belgium
T +32 56 629 531
E info@ckpl.eu

WWW.CKPL.EU

CEKA Multi / CEKA Multi O

CEKA Multi Composite (colour-neutral) and **CEKA Multi O Composite** (opaque-coloured) are self-curing fixing composites. They are suited to fixing zirconium constructions on restorations made from zirconium dioxide, precious metals, non-precious metals and ceramics, by means of adhesion.

CEKA Multi Primer is used as a bonding agent between zirconium dioxide, metals, full ceramic restorations and facing materials.

Let the paste reach room temperature before applying it. Lower temperature can make pressing out and mixing difficult. Quickly use the composite after pressing it out from the automix tip.

Store CEKA Multi/Multi O Composite at 3-9 °C (37-48 °F).
Store CEKA Multi Primer at 10-25 °C (50-77 °F).

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