



SAM CHEST SEAL TRAINING

SAM[®]
MEDICAL

SAM[®]
MEDICAL

WHY WE NEED CHEST SEALS

WHY WE NEED CHEST SEALS



Open chest wounds, specifically Tension Pneumothorax are the leading cause of potentially survivable battlefield deaths.

WHAT IS AN OPEN CHEST WOUND?

WHY WE NEED CHEST SEALS



Open chest wounds, specifically Tension Pneumothorax, are the leading cause of potentially survivable battlefield deaths. A valved chest seal should cover and protect the open chest wound while not preventing potentially lethal pressure from escaping the pleural space.

1. Eastridge BJ, Mabry RL, Seguin P, Cantrell J, Tops T, Uribe P, Mallett O, Zubko T, Oetjen-Gerdes L, Rasmussen TE, Butler FK, Kotwal RS, Holcomb JB, Wade C, Champion H, Lawnick M, Moores L, Blackbourne LH. Death on the battlefield (2001-2011): implications for the future of combat casualty care. *J Trauma Acute Care Surg.* 2012 Dec;73(6 Suppl 5):S431-7.

2. Holcomb JB, McMullin NR, Pearse L, et al. Causes of Death in U.S. Special Operations Forces in the Global War on Terrorism: 2001-2004. *Annals of Surgery.* 2007;245(6):986-991.

A SWAT officer in tactical gear is shown applying a SAM Chest Seal to a person lying on the ground. The officer's uniform has a "SWAT" patch. The seal is a circular, adhesive device with a central white cap and a mesh-like surface. The scene is dimly lit, suggesting an indoor or nighttime setting. In the background, a piece of equipment is labeled "E-100 PERFORMANCE FOC GENERATOR".

SAM CHEST SEAL

<https://www.youtube.com/watch?v=H4RRmzjrsB0>

A man is lying in a hospital bed, appearing to be in pain or unconscious. He has a SAM medical device (a large, adhesive, multi-lead electrode pad) attached to his chest. The device is covered in a red, gel-like substance. A person's hands are visible, one resting on the man's shoulder and another near the device. The background shows the metal frame of the hospital bed. The overall scene is dimly lit, emphasizing the medical emergency.

ENGINEERED TO PRESERVE LIFE

ENGINEERED TO PRESERVE LIFE



ENGINEERED TO PRESERVE LIFE



TRUFLOW™ VALVE*

Is composed of a reinforced dome which protects the chest wound and the internal one-way valve from external forces allowing air to flow from the chest cavity and preventing air ingress into the chest cavity.

*Does not come with non-valved version

https://www.youtube.com/watch?time_continue=2&v=IV17aBtQslg

ENGINEERED TO PRESERVE LIFE

FORTIFIED GEL

Strong Hydrogel tightly adheres to the chest wall in the presence of blood, hair, sweat or sand. And if necessary allows lifting and re-application while maintaining superior adherence.



REVERSE SIDE

ENGINEERED TO PRESERVE LIFE

THIN EDGE LAYER

The thin edge layer resists lifting or peeling, maintaining a formidable seal



ENGINEERED TO PRESERVE LIFE

NIGHT VISION OPTIMIZED

Optimized for night vision, the packaging, pouch, dressing backing, and tabs provide ideal contrast for visibility in normal light, low light, and night conditions (using a night vision device).



ENGINEERED TO PRESERVE LIFE

OVAL SEAL

Larger than other dressings, its oval shape maximizes the surface area of the seal, and its transparency, elasticity and conformability ensure effective adhesion.



ENGINEERED TO PRESERVE LIFE



DUAL APPLICATION TABS

Tabs facilitate the proper placement and removal of adhesive seal, with the ability to vent and reseal.

ENGINEERED TO PRESERVE LIFE

PEEL / TEAR PACKAGING

Quickly and easily open SAM[®] Chest Seal packaging by peeling from its tabs or tearing along the chevron-marked areas.



ENGINEERED TO PRESERVE LIFE



ELASTICITY

Flexible polyurethane film allows the SAM Chest Seal to stretch and contract with normal respiratory movements, greatly improving adherence to the chest wall.

FOUR CONFIGURATIONS

FOUR CONFIGURATIONS



COMBO



VALVE 2.0



NON-VALVED



WITH VALVE

COMBO



NON-VALVED



VALVED 2.0



WITH VALVE






WHY SAM CHEST SEAL?

SAM CHEST SEAL

SAM[®] Chest Seal is engineered to treat, seal, and reseal open chest wounds under most circumstances. After simple and quick application, SAM[®] Chest Seal sticks relentlessly – through extreme heat or cold, no matter the elements. Meets CoTCCC-preferred features.

WHY SAM CHEST SEAL?



**OCCCLUSION -RESISTANT
TRUFLOW VALVE**

WHY SAM CHEST SEAL?

STERILE R

STERILE SEAL


Minimize further
contamination of the wound

WHY SAM CHEST SEAL?

FORTIFIED GEL

Strong Hydrogel adheres in the presence of blood, hair, sweat and sand, or when submersed in water.

WHY SAM CHEST SEAL?



LARGE OVAL SIZE

Larger than other dressings, its oval shape maximizes the surface area of the seal, and its transparency, elasticity and conformability ensure effective adhesion.

WHY SAM CHEST SEAL?



NIGHT VISION OPTIMIZED

Optimized for night vision, the packaging, pouch, dressing backing, and tabs provide ideal contrast for visibility in normal light, low light, and night conditions (using a night vision device).

WHY SAM CHEST SEAL?



6-YEAR SHELF LIFE

Industry-leading shelf life ensures SAM[®] Chest Seal is always ready when you need it.

WHY SAM CHEST SEAL?



A person in military uniform is lying on a patterned rug. A SAM chest seal is being applied to their chest. The seal is a rectangular piece of fabric with a central opening and a strap. The person's hands are visible, holding the seal in place. The background is dark and out of focus.

SAM CHEST SEAL APPLICATION

COMBO



https://www.youtube.com/watch?v=y1log_cklmM




STEP 1

Quickly and easily open the sterilized SAM[®] Chest Seal packaging by peeling from it tabs or tearing along the chevron-marked areas.




STEP 2

Clean the wounded area.




STEP 3

Next, grip the application tab and remove the clear liner.




STEP 4

Center the TRUFLOW Valve over the wound and place the dressing adhesive side down.



STEP 5

Press firmly across its entire surface to ensure adhesion.




STEP 1

**Locate the wound and
remove the casualties
clothes.**



STEP 2

Clean the wounded area.





STEP 3

Next, grip the application tab and remove the clear liner.



STEP 4

Center the dressing over the wound and place the adhesive side down before pressing firmly across its entire surface to ensure adhesion.



A medical professional is performing a procedure on a patient's chest. The patient is lying on a gurney, and the professional is using a SAM (Seal, Apply, Mangle) device to seal a wound. The SAM device is a white, circular, adhesive dressing that is being applied to the patient's chest. The professional's hands are visible, and they are carefully positioning the SAM device over the wound. The patient's chest is the focus of the procedure, and the SAM device is being used to seal the wound. The professional is wearing a white coat, and the patient is wearing a white shirt. The background is a hospital setting, and the overall scene is dimly lit.

SUMMARY

SUMMARY

STERILE

SUPERIOR ADHESION

OCCLUSION-RESISTANT TRUFLOW VALVE

LARGE SIZE

NIGHT VISION OPTIMIZED

6 YEAR SHELF LIFE



FAQS

FAQS



Why is the SAM Chest Seal a large oval shape?

The SAM Chest Seal is oval shaped and larger than other dressings. This maximizes the surface area of the Seal to assure effectiveness of adhesion. If the seal displaces from chest, it cannot protect the wound or the the air flow from chest. The large oval shape also facilitates coverage of the chest wound and surrounding tissue.

FAQS

Why does the SAM Chest Seal come with and without a valve?

The SAM Chest Seal comes in both versions because not all chest injuries are treated the same. The “valved” SAM Chest Seal provides a one-way valve to allow air flow out of the chest thereby preventing the life threatening condition; tension pneumothorax.

The non-valved version provides a large occlusive dressing. Primarily used for chest wounds when medical professionals have the means to carefully monitor for and treat tension pneumothorax using more definitive means.





FAQS

Can the SAM Chest Seal be used in the prehospital setting and austere conditions?

Yes, the SAM Chest Seal has been designed to work in many different environments including extreme temperatures (hot or cold), wet (blood, sweat, rain, etc.), in the presence of body hair, and dirt (sand, dust). The SAM Chest Seal has been tested to assure you, the rescuer that even in truly austere conditions the SAM Chest Seal will work effectively and efficiently each and every time.

FAQS



Why are there so many notches in the SAM Chest Seal packaging?

The notches at the four corners allow quick access to the SAM Chest Seal in less than optimal conditions. No more searching for the tab or being forced to cut the package open. The top flap opening allows for the SAM Chest Seal to be opened such that the packaging itself provides a sterile field. In the field or in the hospital, the packaging of the SAM Chest Seal allows the most effective access when seconds count.

REGULATORY

EC REP


CE 0050



Rx ONLY

STERILE R

Emergo Europe, Prinsessegracht 20,
2514 AP, The Hague, The Netherlands



SAM Medical Products®
27350 SW 95th Ave, Ste 3038
Wilsonville, OR 97070 USA

MORE INFO

sammedical.com

THANK YOU

SAM[®]
MEDICAL

MORE THAN SURVIVAL™