**Uriflex®**  
is a brand of  
**Urifon BV**

Mennigtenbuurt 29  
1427 AX - Amstelhoek (NL)  
Tel: +31 0 2977 - 760 001  
Fax: +31 0 2977 - 760 002

**ApollO24**  
DISCRETE DURING THE DAY  
LOUD AT NIGHT!

Connecting cord orange clip (Dry-Mate inlay) /  
Anschlusskabel orange Clip (Dry-Mate Einlage) /  
Cordon de connexion clip orange (protection Dry-Mate) /  
Cavo di collegamento clip arancione (Inserito Dry-Mate) /  
Cable de conexión clip naranja (compresa Dry-Mate) /  
Przewód z pomarańczowym zaciskiem (wkładki Dry-Mate)

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**CONGRATULATIONS ON YOUR PURCHASE OF THE APOLLO 24!**

The Apollo 24 has many benefits: It is small, light and features a vibration and variable setting audible alarm. At the same time you will receive 2 free inlays which we recommend you use (with the orange clip). This gives you a high sensor sensitivity and the urination episode will be discovered sooner. The alarm will go off sooner and thus enabling training to proceed more efficiently. More inlays can be ordered at low cost from your sales outlet. Before using the alarm, you will need to fit the battery first! See "Getting ready for use". We wish you every success with your training! The Uriflex team.

**How do you intend to use the Apollo 24:**  
1. Combined with your own underwear! Use the cord with the blue clip!  
2. We recommend you try the 2 free Dry-Mate disposable inlays included. Use the cord with the orange clip! The high sensitivity of the sensor means that the urination episode will be highlighted sooner. Due to the sharp teeth in the orange clip, do not attach it to your own clothing!

All components are also available separately at [www.uriflex.com](http://www.uriflex.com) or from your local sales outlet.

**Getting ready for use**  
Release the 2 battery cover screws using the screwdriver included. Fit the battery with the text upwards (see illustration). Refit the cover and tighten the screws. The red indicator light on the front will flash briefly once every 20-30 seconds.

**Setting alarm signal type**  
You have a choice of vibration, sound or both combined. The standard setting for the alarm is sound only. If you wish to change this, press and hold the Stop button (see "illustration") for 3 seconds, until the red indicator light is on constantly. Press button A or B. Each time you press, the setting will change from sound only alarm, vibration only or combined sound and vibration. Press the stop button once you have reached the desired alarm type. The red indicator light on constantly will go out. Note: not the red indicator light must be on constantly if you wish to set the alarm type!

**Setting the desired sound signal volume**  
There are 3 settings: soft, medium and loud. The factory setting is medium. Pressing button A for at least 1 second increases the volume. Pressing button B lowers the volume in steps.

**Low battery signal**  
If the red indicator light is on continuously of its own accord, it means the battery is almost flat. Replace the battery.

**Using the Apollo 24**  
Attach the clip over the pyjama collar of fit the Apollo 24 with a safety pin wherever you like. Pass the cord under the clothing to the underpants. Press the Stop button when the alarm goes off in order to switch it off. The alarm will go off again after 20 seconds if the connecting cord is still attached at both ends. This will continue until the connecting cord is detached at either end.

**Then the alarm goes off...**  
Urinating causes the alarm to go off. The user will wake up and the bladder sphincter will close by reflex action. The user gets out of bed, detaches the connecting cord and switches the alarm off with the Stop button. Urination is completed in the toilet. The cord is then reattached to a new Dry-Mate inlay (with an orange clip) or directly to the underpants (with a blue clip) and the user goes back to sleep.

**Important: All in full awareness!**  
It is important that the user should be well awake and fully aware of everything. It may be desirable for a parent, guardian or carer to sleep in the same room as the user, in order to ensure that there is an immediate reaction to the alarm.

An immediate reaction to the alarm signal is the key to success. Only then will the reflex action of "timed" waking be encouraged and developed further. That is why we recommend the use of Dry-Mate inlays. The very high sensitivity of the sensor will react to each drop of urine, whereas the small clip on its own often requires more urine to activate the alarm.

By using the Apollo 24 on a daily basis, after a couple of weeks the user will wake up of his or her own accord and in time.

**Fill in the score card included every morning**  
Every morning the score card together with show how many points were earned that night. This will naturally provide a great deal of encouragement. This will often work even better, if you attach your own - small - reward to it.

**Training successful**  
Stop after 3 weeks of dry nights (no urination at all or woken up of one's own accord) with the bedwetting alarm. If the bed is still dry two weeks after that, then the training has succeeded.

**Use the Apollo 24 as follows.**  
1. **HIGHLY RECOMMENDED: USING APOLLO 24 WITH DRY-MATE INLAY**  
The orange clip to be used on the top edge of the inlay. It is important to place the metal teeth on the clip precisely where the black lines start! The inlay is placed in the user's own underpants using the adhesive edge.

**Benefits**  
+ High sensor sensitivity hence rapid alarm response to a urination episode.  
+ To be used in user's own underwear.  
+ Does not leak, hence no need to launder underwear if reacted to quickly.  
**NB!** The orange clip to be used only in combination with the Dry-Mate inlay on account of sharp teeth! Not to be used directly on user's own clothing. Use the blue clip for this!

**2. USING APOLLO 24 WITH USER'S OWN UNDERWEAR**  
The blue clip is to be attached to the outside of the user's own underpants, at the point where you expect the first drops of urine.  
Note that the fabric in the clip must be wet in order to activate the alarm, so consider carefully where you place it. The sooner the urine loss is spotted, the more effective the training will be.

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