

PRODUCT DATASHEET

C3 Emptying Bag w/SAP by Curas™	
Brand	: CURAS®
Category	: Urology
Range	: Drainage Bags
Product REF	: 20465
Product Description	: C3 Emptying Bag w/SAP by Curas™, 1.5L, Dual-NRV System (incl. Spigot), Medium Flexi-Funnel™, Cross-Hanger™, Non-DEHP, CleanPak™, WardPak™



PRODUCT DESCRIPTION

The C3 Emptying Bag by Curas incorporates our unique innovation of a clean disposable and exclusively closed urine drainage system. The C3 Emptying Bag is designed for a safer emptying and drainage of the patient's urine, when removed from the catheterized patients' urine drainage bag.

More importantly, C3 Emptying Bags come with the Medium Flexi-Funnel™ feature that allows safe and easy connection to all forms of bottom outlets/drainage taps or tubes, irrespective of whether the patients' drainable urine bag is a regular urine bag or a more advanced urine meter system.

This C3 Emptying Bag comes included with SAP by Curas™ (Super Absorbent Polymer) for immediate solidification and safe handling of urine. The C3 Emptying Bags are also fitted with Dual-NRV System (including Spigot for SAP models), as a default solution to further avoid risks of spillage. The product comes with a 1.5L capacity.

REF	Description	Capacity
20465	C3 Emptying Bag w/SAP by Curas™, Dual-NRV System (incl. Spigot), Medium Flexi-Funnel™, Cross-Hanger™, Non-DEHP, CleanPak™, WardPak™	1.5L

INTENDED USE

As depicted by its name, the unique Curas C3 Emptying Bags are used for draining and emptying urine from patients' drainable bags and/or hourly urine measurement systems (urine meters). Designed to assist nurses during critical and highly infectious conditions. It is a much safer, cleaner, convenient, and a hygienic option in comparison to the commonly used jug, urinal, or otherwise open emptying method.

Diagram below illustrates how the Curas C3 Emptying Bag should be connected to a urine meter.



FEATURES & BENEFITS

- The C3 Emptying Bag provides a safer and more convenient method of draining and emptying urine from the patient's drainable urine bag.
- Super Absorbent Polymer (SAP by Curas™) rapidly crystallizes urine; thus, ensuring safe disposal of urine without any risk of spillage.
- Clear ml-graduation markings on C3 Emptying Bag allow easy reading on the approximate amount of urine collected.
- The Medium Flexi-Funnel™ fits all standard bags (e.g. C2 bag) and is compatible with Urine Meters (e.g. B. Braun-Ureofix 500 & Bard urine meter drainage bag).
- The Cross Hanger™ is a uniquely designed cord system that provides extra security for the nurses when hanging, holding and carrying the C3.
- Dual NRV System (including spigot) ensures no reflux of urine and risk of urine spillage thus ensuring the safe transfer of C3 bag for disposal.
- Double-sealed bag and premium level foil material are specially designed to provide extra strength during use.
- CleanPak™ is designed to provide 100% premium protection of hygienic integrity.
- C3 packaging comes with data (LOT, EXP), and product tracking via bar-coding for added control.
- Medium-level packaging system, WardPak™ dispensing box, allows easy use at ward level.

PRODUCT SPECIFICATIONS

REF	SAP	Capacity	Non-Return Valve	Spigot	Flexi-Funnel™	Cross Hanger™	CleanPak™ (Individual Inner Bag)	WardPak™ (Box)	Weight/pc
20465	Yes	1.5L	Yes	Yes	Yes (Medium)	Yes	Yes	Yes	~42g

Composition

- Bag and connector made of non-DEHP based medical grade PVC.
- Cord made of Dacron.
- Spigot made of PE.
- Super Absorbent Polymer (SAP by Curas™) composed of sodium polyacrylate powder.
- All components are latex-free.
- Material contains di(2-ethylhexyl) phthalate (DEHP) (Bag:12ppm)

Country of Origin

- China

Sterilization

- Non-sterile

Shelf Life

- 5 years

Selected Standards & Guidelines

- | | | |
|-------------------------|--------------------|-----------------------|
| ■ ISO 9001:2015 | ■ ISO 14971:2019 | ■ EN ISO 14698-1:2003 |
| ■ EN ISO 13485:2016 | ■ ISO 20417:2021 | ■ EN ISO 14698-2:2003 |
| ■ MDR 2017/745, Class I | ■ ISO 15223-1:2016 | ■ EN ISO 10993-1:2018 |

Recommended Storage Condition

- Store in dry, clean, and well-ventilated area at room temperatures between 15°C to 25°C (59°-77°F).
- Avoid exposure to extreme temperatures. Keep away from moisture, heat and sources of ignition.
- Do not expose the product to radiation, sunlight and UV light, as it will change the properties of the product.

Proposed Disposal Method

- Dispose of used product into clinical waste bag or adhere to the clinical waste protocol, e.g., incineration

Product Symbols







Symbols for SAP by Curas™ only



PACKAGING INFO

REF	GTIN	CleanPak™	WardPak™	Carton	Gross Weight
20465	05051400204654	1 pc / bag	25 pcs / box	200 pcs / carton	~13 kg

Components / Material Properties			Recommendations on Waste Handling			
Packaging	Material	Symbol	Incinerate	Landfill	Recycle	Reuse
CleanPak™	HDPE		√	(√)		
WardPak™	Corrugated Board		(√)	(√)	√	
Outer Carton	Corrugated Board		(√)	(√)	√	
Transport Pallet	Wood				(√)	√

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- HS Code: 3926909790
- Documentation available upon request. Kindly contact: info@curas.com