
optiflushTM SP

Pre-filled Flush Syringe

Instructions for use

GBR

OptiFlush SP Pre-filled Flush Syringe

PRODUCT REF / CONTENTS:

Product Code / Description

5003	3ml syringe
5005	5ml syringe
5010	10ml syringe

PRODUCT STRUCTURE: (See figure 1.)

1. Protective Tip Cap - Luer Lock
2. 0.9% sodium chloride
3. Syringe Barrel
4. Piston
5. Plunger Rod

INTENDED USE:

The device is used for closing and flushing the end of catheters at intervals of different drug treatment.

CONTRAINDICATIONS: None.

PRECAUTIONS:

1. The use of the product must meet the requirements of relevant operational specifications and relevant regulations of medical department, and be restricted to trained doctors or nurses.
2. Check for package integrity. It is strictly prohibited to use if product package is damaged, or the syringe breaks, leaks, disintegrates.
3. It is strictly prohibited to use if protective tip cap is not installed correctly, loose or falls off.
4. Visual examination for clarity of solution, if solution discolouration, turbidity, precipitation or any form of suspended particulate matter material appears-DO NOT USE.
5. DO NOT RE-STERILIZE
6. Please use the product within the validity period. The product batch number, production date and expiration date are shown on the product packaging label.
7. It is strictly prohibited to re-use. Single use only, and discard any unused remaining solution.
8. Check the compatibility of 0.9% sodium chloride solution with drugs. Direct contact with drugs incompatible with 0.9% sodium chloride is strictly prohibited.
9. It is strictly prohibited to be used for subcutaneous and intramuscular injection and dispensing. Only for flushing of intra-vascular catheters. It is not for medication dilutions, nor for dry product reconstitution.
10. DO NOT PLACE THE PRODUCT ON A STERILE FIELD.
11. Do not use if you are allergic to the sodium chloride solution.

Any serious incident or malfunction that has occurred in relation to **OptiFlush SP** should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

SHELF LIFE: 3 Years

PRODUCTION DATE: Refer to the package label.

STERILISATION METHOD: Steam sterilisation.

STORAGE: Store at controlled temperature (10-30°). Do not freeze.

DIRECTIONS FOR USE: (See figure 2.)

(FOLLOW ASEPTIC TECHNIQUE THROUGHOUT THE PROCEDURE)

1. Open syringe package and remove syringe, inspect clarity of solution, check the protective tip cap is in place.
2. Push the plunger rod up to release the seal between the piston/plunger and the barrel before removing the protective tip cap.
3. Remove protective tip cap by twisting off aseptically.
4. Expel air by pushing plunger forward.
5. Connect flush syringe to catheter aseptically and push the plunger rod to flush the required volume according to the institution's policy.
6. Discard the used syringe and residual part of solution. DO NOT RE-USE.

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


Jiangsu Embrace Science&Technology Development Co., Ltd.
29-2 Yingbin Road, Huai'an Economic and Technological
Development Zone, Jiangsu, P.R.C.

Sponsor:
Optimum Medical Australasia Pty Limited,
Level 4, 470 Church Street,
North Parramatta New South Wales, 2151, Australia

EU REP

Riomavix S.L.
Calle de Almansa 55, 1D,
Madrid 28039, Spain

UKRP

Rociale Healthcare Limited
Cwm Cynon Business Park,
Mountain Ash, Rhondda,
United Kingdom, CF45 4ER

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 - Suite 123,
2595 AM, The Hague, The Netherlands

	Legal manufacturer		Importer - Indicates the entity importing the medical device into the European market
	Authorised Representative in the European Community		Distributor - Indicates the entity distributing the medical device into the locale
	Indicates that the product is in compliance with European legislation for medical devices		Medical Device
	UDI - Indicates Unique Device Information		Use-by Date YYYY-MM-DD
	Date of Manufacture		Catalogue Number
	Batch Code		Caution
	Country of manufacture - China		Consult Instructions for Use
	Keep dry		Keep Away from Sunlight
	Latex Free		Do not Use if Damaged and Consult IFU
	Single use on a single patient		Do not resterilise
	Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile.		Indicates the upper and lower limit of temperature to which the medical device can be safely exposed
	Demonstrates the packaging layer plastic type (resin code 5/PP)		Demonstrates the packaging layer plastic type (resin code 4/LDPE)
	Demonstrates that the packaging layer is capable of being recycled		Demonstrates that the paper elements of the packaging only are capable of being recycled

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Optimum Medical Solutions Ltd