Digital Microscope with LED, Manual

SA-9757





Specification

Image CMOS Sensor Controller High Speed DSP (Driver Free Available) Focus range 10 mm ~ far sight Snap Shot Software and Hardware Video capture resolution Standard 640*480, Max 1600*1200 Built-in 8 white-light LED and adjustable Still image resolution 640*480, Max 1600*1200 Digital Zoom 5X Brightness control Power supply USB port (5V DC) Dimension 112 mm (L) X 33 mm (R) Net weight of product: 250 g

Index

- 1. Specifications
- 2. Table of contents
- Package Contents Tips and warnings Features and Components
- 4. Software installation PortableCapture-Plus Windows
- 5. Software installation PortableCapture-Plus MAC OS
- Setting software Language Menu options
- 7. Data speed Photo Capture
- 8. Video recording
- 9. Timer
- 10. Image editing
- 11. Photo editing Drawing
- 12. Text Measurements
- 13. Calibration
- 16. Adjust calibration

Icons and statements 17.

Content

- 1. Digital Microscope
- 2. USB 2.0 cable
- 3. Metal fixed stand
- 4. CD ROM (driver, measurement software, user's manual)

Tips and warnings

- 1. Do not use the LCD digital microscope in environments where such devices are prohibited, improper use may result in serious accidents.
- 2. Protect the LCD digital microscope from water and other liquids. Do not use the LCD digital microscope with wet hands. Never use the LCD digital microscope in the rain or snow. Moisture can cause fire and electric shock. If you notice smoke or odors coming from your LCD digital microscope, turn it off immediately.
- 3. The LCD digital microscope may become slightly warm during use. This is not a malfunction.
- 4. Wipe the LCD digital microscope with a soft, dry cloth.

Function - Parts



- 1. USB Microscope
- 2. Focus
- 3. LED
- 4. Lens
- 5. Connection cable USB
- 6. Brightness control
- 7. USB plug

Software installation PortableCapture-Plus Windows

Place the CD in the CD-ROM drive.

The Capture-Plus software starts up automatically for Windows users. If the installation does not start automatically, select "setup.exe" on the CD, to install the software on the computer.

Organisieren • Auf Datentra	ger prennen		822 •
Name	Änderungsdatum	Тур	Größe
Momentan auf dem Daten	träger vorhandene Dateien (7)	
🎍 mac	19.11.2013 08:44	Dateiordner	
Ja src	19.11.2013 08:44	Dateiordner	
퉬 User's manual	19.11.2013 08:44	Dateiordner	
AutoRun	14.11.2013 08:20	Setup-Informatio	1 KB
Menu Menu	14.11.2013 08:20	Anwendung	359 KB
🜉 Setup	14.11.2013 08:20	Bitmap-Bild	4 KB
🕝 setup	14.11.2013 08:20	Anwendung	9.981 KB
Dateien, die auf den Daten	nträger geschrieben werden s	sollen (1)	
a desktop	04.05.2015 15:07	Konfigurationsein	1 KB



Then click Install MicroCapture, the installation is completely automatic.

Reset the computer after installing the driver.

Installation may be blocked by the antivirus software, so have the antivirus software allow all processes or close the antivirus before installing.

If the installation was successful, you will see the icon on the toolbar, double-click it to run it.



Software installation PortableCapture-Plus MAC OS

Place the CD in the CD-ROM drive.

For MAC OS, there is a MAC folder on the CD.

000	i mac			C
e	1 Objekt, 0 KB verfügbar			
Name	Anderungsdatum	Größe	Art	
MicroCapture_PLUS_Package.pkg	14. November 2013 08:20	20,9 MB	Installspaket	

Copy the MAC file MicroCapture-Plus to the desktop. Double-click on the copied file to open it. The installation wizard will launch and guide you through the rest of the process.

For faster and easier access to the software, we recommend that you create a shortcut to the program after installation. After installation, we recommend that you create a shortcut to the program and place it on your desktop.

Reset the computer after installing the driver.

Installation may be blocked by the antivirus software, so have the antivirus software allow all processes or close the antivirus before installing.

If the installation was successful, you will see the icon on the toolbar, double-click it to run it.



Setting up software

Launch the software by clicking on the icon



The following menu will appear on your screen



Language

By default, the software is set to English. To change the language, select Options > Language > Dutch.



Menu options

Open the Options menu and the following screen appears



- 1. Resolution Selection of image resolution
- 2. Date/Time Display/hide date and time,(is not saved on the photo/video)
- 3. Language Selection of the menu language in the software
- 4. Show/hide reticle

Data speed

The transfer rate on the USB interface depends on which version you have 2.0 or 3.0 Therefore, for smooth image display, the resolution is preset to 1280 x 960. Of course, photos and videos are saved at the resolution selected under Options > Resolution.

Photo Capture

To take photos, choose from the menu below



Then click on Photo



After this, you can press the record button located in the connection cable. The recording is saved in the File > Photo folder in the main menu. At the same time, the photos are displayed as thumbnails on the right side of the program window



When you click the mouse on a thumbnail, the selected photo is displayed in the larger view. If you right-click on a thumbnail, you will see a submenu, which allows you to perform different actions for the selected photo.



- 1. Save
- 2. Сору
- 3. Clear each
- 4. Clear all

Video recording

To create a video, choose from the menu below



Then click on Video



The above icon starts flashing red, to stop recording click the video icon again and the flashing will stop.

The recorded videos are saved on your computer in the Microcapture_Video folder. You can also access this folder via File > Video folder in the main menu.



When you click the mouse on a thumbnail, the selected video is displayed in the larger view. If you rightclick on a thumbnail, you will see a submenu, which allows you to perform different actions for the selected video.



- 1. Save
- 2. Copy
- 3. Clear each
- 4. Clear all

Timer

You have the option of setting a timer to record photos and videos. To do this, click on the icon



the following window opens

🔬 Zeitgeber	×		
Foto	Video		
Startzeit	2015/10/21 10:50:01 -		
Interval(Sek.):	10		
Videoaufnahmezei	t(Sek.): 10		
Anzahl:	10		
Erfasste Anzahl:			
Start	Stop		

Photo/Video - Select whether to record photos or videos. Start time - Set recording start time Interval (sec.) Set recording interval (from 0 to 99 seconds) Video recording time (sec.) Set length of video recording (from 0 to 99 seconds) Number - Set number of shots (from 0 to 99)

You can see how many recordings have already been made in the Recorded Number field.

Image editing

The software also has an image editing program. To edit a captured photo, open it by double-clicking on the thumbnail or right-clicking on the thumbnail and then selecting Open. The photo will open in a separate window.



Photo editing

You now have several options for editing the photo



Save photo



Copy photo to clipboard



Copy photo



Previous



Next



Undo



Again

Drawing

Using the icon below, you can select drawing tools under the pull-down menu as below.



Selection options for drawing; Line, Arrow, Free hand, Rectangle



Using the icon below, you can select font under the pull-down menu as below.



Selection options for text input; Font, Text color



Measurements

Using the icon below, you can select measurement under the pull-down menu as below.



Selection options for text input; Angle (simple line with one start and one end point), Continuous line, Radius, Diameter, 3 Point angle, Line properties, Text properties, Unit of measure (mm, cm, pixel, inch).



The following setting/measurements are from Left to Right;

Line, ring, Square, Polygon, Line length, Oval, Arc, Point to point distance, Point to line distance, Parallel line, Angle.

1. Line measurement at any angle Simply click the mouse button on one point and drag to the next point, then release.

2. Continuous Line Measurement With the continuous line measurement, you can measure any continuous length you want. Simply click the mouse button on one point and drag to the other desired point and then release.

3. Circle radius measurement Simply click and extend to the desired radius, the radius, circumference and area of the circle will appear accordingly.

4. Circle diameter measurement Simply click and extend to the desired diameter, the diameter, circumference and area of the circle will appear accordingly.

5. Three-point angle measurement Just define any three points and click, the value of the angle will then appear.

Calibration

The software provides a calibration function that allows you to perform the most accurate measurements on your photos. You should always set this function first if you want to obtain correct value when taking measurements.



You can use the Crossing option to check that the scale unit is parallel. Then adjust the magnification until the image is clear. Take a picture of the selected area when the scale unit is parallel to the crosswire, as below.

Our experience taught us that the measurement could be controlled with 99% high precision. Measurement tolerance can be caused by distortion and blur of the human eye.



Double-click the image to open the measurement interface, as below.



Then click on the calibration icon



Click on the start and end points over a known value. Now a dialog box like below will appear and you need to enter the exact value of the trial measurement in Actual measurements.

For correct measurement use the enclosed Calibration Ruler



Now manually enter the actual size in mm of the distance you just defined in the field.



The software will then automatically calibrate the magnification. Click OK and the adjusted magnification factor will appear in the upper right corner. Reliable measurements can now be made on all photos taken after calibration using the same microscope focus setting. If you change the focus for

a photograph using the dial on the device, it is necessary to recalibrate. If necessary, you can also manually set the magnification for a single photo or for future shots.

You can now use any line option from the menu under the icon below



Use to measure the value again. If the value measured with software is close to or exactly 100% equal to the exact value of the calibration ruler, it means the calibration is complete.



Adjust calibration

To change the magnification and adjust the focus for photo shooting, recalibrate the magnification factor.



Reset magnification > Reset picture magnification. And then restart calibration (see earlier in this manual)



Change the magnification factor and take several pictures, click on the image to open the measurement interface.

And click on the icon below



Set magnification > Set picture magnification, now a dialog box will appear, as below;



Enter the magnification factor in the dialog box and click OK.

If you want to use the same magnification factor to take pictures, just click on



Set magnification > Set capture magnification.

Now all images will use the same magnification factor and you can start measuring immediately.



Icons and the explanations



Take picture - Used to take a picture one at a time.



Start/Stop recording video.



Timer for photo/video - Set the start time, interval time and total number of photos/videos.



Full screen mode Press "ESC" or double-click to exit.



Save As



Copy to clipboard



Copy file



Previous



Next



Undo



Again



Drawing



Text



Measurement



Calibration



View, pause or stop recording.

Your supplier

