

drive

 **DeVilbiss**[®]
HEALTHCARE



X-Fold Rollator

Instruction Manual

Please read these instructions carefully before assembling or using this product.

CE

product description

Intended User Group

The intended user of the X-Fold Rollator is any person who has difficulty in walking and / or who has reduced mobility.

Intended Purpose

A walking device manipulated with both arms to assist a person to maintain stability and balance whilst walking.

Environment of use

- Indoor use in domestic, residential and / or hospital environments.
- Outdoor use on firm, even surfaces.

Indications

To provide support to a user who has difficulty in walking and / or who has reduced mobility.

Contraindications

- Unable to weight bear
- Unable to adequately manoeuvre the rollator or apply the brakes
- Users with Hemiplegia
- Users with only one arm
- Users who weigh in excess of 135 kg
- The rollator being used as a transport device

Other contraindications may be relevant which are specific to the patient, or care environment.

specification

Length: 68 cm (27") | Width: 64 cm (25")

Height: 81 – 94 cm (32" – 37")

Seat dimensions (cm) : 23 (L) x 46 (W) x 56 (H)

Seat dimensions (in): 9" (L) x 18" (W) x 22" (H)

Mass: 8.35 kg (18.4 lb)

Mass Limit: 135 kg (21 stone)

Bag Mass Limit: 5 kg (11 lb)

symbol definition

The following symbols are found on this product:



Conforms to the Medical Devices Regulation
EU 2017/745 (see declaration of conformity)



Caution - beware
of potential hazard



Manufacturer



Consult
instructions for use



Date of manufacture



Product code



Serial number



Medical device



Quantity

contact information

For assistance in setting up, using, maintaining your rollator, to report unexpected operation or for any service, warranty, sales or customer service information regarding this product, please contact your provider.

If in doubt, contact Drive DeVilbiss Healthcare Ltd. using the contact details on the back page.

If use of this device results in a serious incident, reports should be forwarded to the manufacturer and local competent authority.

Please quote the product serial code on all correspondence.

Service & Maintenance

Customer Service

Tel: +44 (0)1422 233136

Tel: +44 (0) 845 0600 333

Fax: +44 (0)1422 233010

Fax: +44 (0) 845 0600 334

Spares

Tel: +44 (0)1274 475125

general warnings

Please read all of these instructions fully before assembly or use.

- The rollator is a walking aid only and must not be used as a transport device (i.e. wheelchair).
- Consult a healthcare professional before using any assisted mobility device. A risk assessment must be carried out to determine the suitability of the rollator for the specific user.
- Do not exceed the 135 kg (21 st) user mass limit.
- Dispose of plastic bags and all other packing, to avoid risk of suffocation to children or animals.
- Make sure the rollator is fully open, the seat is down and the frame is locked securely before use.
- Ensure that hand grips are at equal heights. Test the rollator for stability before applying full weight to the hand grips.
- Before every use, check the brakes stop the rear wheels moving. Do not use a rollator with damaged or faulty brakes.
- Do not use the rollator on gradients above 7°.
- Do not use on extremely soft, uneven or wet ground, where the rollator may lose traction.
- Do not use a damaged / badly worn seat.
- Always lock the brakes before sitting down. Do not use the resting seat on inclined surfaces.
- Take care not to trap body parts in gaps around the rollator during use.
- Do not place the rollator near a heating source, as it may cause surface temperatures to increase.
- Do not use the rollator near ignition sources such as cigarettes, matches, fireplaces etc.
- Ensure that good health and safety practice is used for lifting and carrying the rollator.
- Replace the tyres when significantly worn or damaged. Contact Drive DeVilbiss Healthcare Ltd for approved replacement parts.

parts description




setup



1. Install the back rest by sliding over the bracket. Use a screwdriver to tighten the screw.



2. Install the top section of the rollator by pulling upwards on the back rest until fully unfolded.



3.

Lock the top section upright by flipping over the camlock levers.

Camlock lever

4.

Secure the camlock levers by pushing the levers in to the slots shown.



5.

Push down on the cross brace rails until the rails lock into the plastic slots.

6.

To attach the carry bag, hook the bag loops on both ends of the cross brace rails under the back rest (as shown).



using the X-Fold

Once the frame has been unfolded and the cross brace is secure, the X-Fold is ready for use. Steer the rollator by holding both handles.



Do not use a damaged or badly worn rollator (see care and maintenance).

Ensure the camlocks are pushed into the side slots before every use (see step 4 of setup).


Height Adjustment

The height can be adjusted to suit users of different heights. Each height setting is numbered from 1 to 6. The handles should rest at wrist height.

1. Remove any heavy items from the rollator.
2. Push the button on one side (as shown), while raising or lowering the handle. Keep moving the handle until the button 'clicks' to another setting.
3. Keep pushing the button until the handle is at the desired setting.
4. Adjust the height of the other handle using the second button. Make sure the buttons are aligned to the same number.



Height
Adjustment
Button



Brake Operation

To use the parking brakes, push both brake levers down (see “A”) until the brakes are locked. To release the brakes, pull both levers up and release.

To use the brakes to slow the rollator down, pull up on both brake levers (see “B”). Release the levers to take the brakes off.

Brake Adjustment

The brakes may require adjusting over time. To adjust the braking force, use a 10mm spanner to tighten / loosen the nut as required, while holding the metal tip in place. Always retest the brakes before use.



Always lock the brakes when sitting on the resting seat.



Kerb and Obstacle Climbing

The front wheels of the rollator can be lifted to help negotiate kerbs and other obstacles.

Press down on the stepper lever (as shown) to aid raising the front castors.

storage and transport



Take care not to trap body parts in gaps around the rollator while folding.

The X-Fold can be folded to save space for storage and transport. Make sure not to store the rollator in damp conditions for an extended period.

Folding the rollator

Reverse the setup procedure to fold the rollator (as described below).

1. Remove the bag.
2. Fold the rollator up by pulling the strap in the centre of the resting sling (as shown).
3. Release the camlock levers and flip over to release the top section. (See steps 3 and 4 from 'Setup'.)
4. Hold the back rest and slowly lower the top section to fully fold the rollator.



cleaning and disinfection

The use of neat bleach, abrasive products or similar surface cleaners is not recommended, as damage may be caused to the cleaned surfaces.



PPE must be worn during cleaning and disinfection to prevent the risk of infection.

Care should be taken not to leave any cleaning residues on the handgrips or seat which may leave them slippery when used.

It is advisable to remove the carry bag and any other attachments from the rollator.

General Cleaning:

- Use a regular household cleaner and a damp cloth. Start with the cleanest parts of the rollator and systematically move to the dirtiest parts.
- Extra care should be taken around areas where excess dirt or dust may gather. The cloth should be changed during the cleaning process if it becomes soiled.
- Wipe down with a clean cloth, moistened with a mild detergent and diluted in warm water (40°C).
- Rinse with cold, clean water and a clean cloth and allow to dry fully before use.

Disinfection:

- Mop up any fluid with paper towels.
- Clean the rollator using the instructions provided in 'General Cleaning'.
- Wipe all surfaces of the rollator with a disinfectant (e.g. 70% IPA solution).
- Allow contact time and wipe with a clean cloth or rinse, in line with the manufacturer's instructions for the disinfectant.
- Allow to dry fully before use.

care and maintenance



Perform daily inspections before use.
Do not use if any damage or fault is found.

- Store the rollator in a dry place and ensure that it is clean and dry before storing.
- Ensure the walker is kept clean in line with the instructions provided in this instruction manual.
- Check moving parts (e.g. wheels and height adjustment) regularly for obstruction. All moving parts should move freely.
- Ensure the brakes are functioning correctly. When the levers are squeezed or locked, the two rear wheels should not be able to turn. If the brakes loosen over time, tighten the brake adjustment mechanism and retest the brake levers before use.
- All fasteners (screws, nuts and bolts) should be securely tightened and checked regularly.
- Check the front castor wheels swivel freely. Withdraw the walker from use if it does not.
- Periodically spray the exposed brake cable with a water dispersing lubricant to ensure that the brakes continue to work smoothly.

Drive DeVilbiss Healthcare Ltd recommends that the rollator is serviced annually by a suitably qualified person. Contact Drive DeVilbiss Healthcare or your local supplier for approved spare parts.

disposal

When the rollator, or any associated packaging reaches the end of its useful life, follow local recycling policies for disposal.

When unpacked for the first time, the cardboard boxes used for the packaging can be recycled at recycling centres.

The X-Fold Rollator consists of aluminium, steel and plastic components. Individual parts can be separated and disposed according to the material.

warranty details

There is a comprehensive twelve-month warranty from the date on which your new walker is delivered. The warranty covers the walker for repairs or replacement during this period. For more detail, please see the warranty conditions below:

1. Any work or replacement part installation must be carried out by an authorised dealer / service agent.
2. To apply the warranty should your walker require attention please contact the outlet from which you purchased the walker.
3. Should any part of the walker require repair or full or part replacement, as a result of a manufacturing or material defect within the warranty period, parts will be supplied free of charge. Note: The guarantee is not transferable.
4. Any repaired or replaced parts will be covered by the balance of the warranty period on the walker.
5. Parts replaced after the original warranty has expired will be covered by a three month warranty.
6. Consumable items supplied will not generally be covered during the normal warranty period unless such items require repair or replacement clearly as a direct result of a manufacturing or material defect. Such items include (among others) upholstery and tyres.
7. The above warranty conditions apply to brand new walkers, whether your walker is covered contact your dealer.
8. Under normal circumstances, no responsibility will be accepted where the walker has required assistance as a direct result of:
 - a. the walker part not having been maintained in accordance with the manufacturer's recommendations
 - b. failure to use the manufacturer's specified parts
 - c. the walker or part having been damaged due to neglect, accident or improper use
 - d. the walker or part having been altered from the manufacturer's specification or repairs having been attempted before the dealer is notified

In the event of your walker requiring attention, please contact your service agent / dealer and give all relevant details so they can act quickly. This warranty is issued free of charge by the manufacturer to the original retail purchaser of the device. This does not affect your statutory rights.

Drive DeVilbiss Healthcare Ltd, Holmfield,
Halifax, West Yorkshire, HX2 9TN, GB

Tel: 0845 0600 333 | Fax: 0845 0600 334

Email: info@drivedevilbiss.co.uk

www.drivedevilbiss.co.uk