

PureLight® Sp02 Sensors

The PureLight Principle

Accurate and dependable SpO₂ readings start with PureLight Sensors



MOMIN NODEN



Red and infrared light sources (LEDs)

are essential to pulse oximetry sensors

- √ The red LED measures deoxygenated hemoglobin
- The infrared LED measures oxygenated hemoglobin

Common red LEDs may produce a secondary spectrum emission (contamination)

- ✓ Contamination impacts the oxygenated hemoglobin measurement
- \checkmark LED contamination shifts the oximeter's calibration curve resulting in inaccurate SpO $_2$ readings below 80%


PureLight red LEDs only produce the pure light spectrum

- ✓ PureLight LEDs do not cause a shift in the oximeter's calibration curve at SpO₂ levels below 80%
- ✓ Nonin PureLight Sensors eliminate variations in readings from patient-to-patient and sensor-to-sensor

Product Highlights

- Reliable Consistent readings from sensor-to-sensor
- ✓ Accurate PureLight LEDs in every sensor
- √ Functional Various options for multiple site applications
- ✓ Versatile Spot checking, short-term and extended monitoring
- Comprehensive Full range of sizes for all patients
- ✓ Comfortable *Trim design, great fit & latex-free*

Oximeter Calibration Curve





PureLight Sensor

Reusable SpO₂ Sensors

Nonin has designed an array of PureLight reusable sensors — available for patients of all sizes (adult, pediatric, & infant) and clinical presentations. The 8000 Series of reusable sensors are ideal for short-term monitoring, stress testing and spot checks — any situation where the risk of cross-contamination is low. A comfortable sensor, each is designed for a specific application and is form-fitted to decrease ambient light interference.



Reusable SpO₂ sensors:

- ✓ Proven PureLight LEDs for accuracy at critical SpO₂ levels
- ✓ Comfortable, self-aligning grips to assure proper placement
- √ Quick and easy application
- **✓** Variety of finger, earlobe and forehead sensors for various monitoring applications

8000R



Reflectance Sensor

Holder

(10 pack w/ 20 adhesive collars)



✓ Cost-effective and durable, backed with our warranty

√ Easy-to-clean for convenience and reliability

*Optional Reusable Sensor Extension Cables are available in 1m, 3m, 6m and 9m lengths

Reusable Flex S with FlexiWra

Nonin's Flex Sensor System combines of a disposable sensor with the econor reusable Flex Sensor, made of a durab matching disposable FlexiWrap adhes ideal for numerous applications included studies. Available in sizes to fit neonatis a great value for patient monitoring FlexiWraps (additional packages of 25)

Reusable Flex SpO₂ senso

- √ Proven PureLight LEDs for accuracy at
- **✓** Durability for repeated use
- ✓ Optimal performance in motion situation

 or property of the property of the performance in motion situation.

 Optimal performance in motion situation
- **✓** Cost-effective convenience
- ✓ Simple cleaning for ease-of-use
- √ Trim shape for comfort and contoured to the contou

Sensor Type	1 meter	2
8000J	✓	
8008J	✓	
8001J	✓	

^{*}Optional Reusable Sensor Extension Cable

Infant F

Adult Flex System





's – Simply Brilliant

pO, Sensors

nics of a reusable sensor. The le silicone, is secured with a ive. This hybrid sensor system is ling extended monitoring and sleep es to adults, the Flex Sensor System . Each sensor is packaged with 25

the convenience and performance

rs with FlexiWrap:

may be purchased separately).

critical SpO₂ levels

neter	3 meter	Extension Cable*		
	✓	✓		
		✓		
		✓		
re available in 1m, 3m, 6m and 9m lengtl				

Disposable SpO, Sensors

Nonin's Flexi-Form® II sensors are ideal for extended monitoring or situations where the risk of cross-contamination is high. The disposable Flexi-Form II sensors feature an inclusive cable, sensor and adhesive wrap to provide patients ultimate comfort, while providing a disposable solution at an attractive price. Each sensor has a simple placement illustration and an attached over-wrap that acts as a stabilizer for securing the sensor in position. These single-use sensors are available for patients of all sizes.

Disposable SpO₂ sensors:

- ✓ Proven PureLight LEDs for accuracy at critical SpO₂ levels
- √ Self-adhesive for extended monitoring and security
- √ Easy to apply for clinicians or patients
- √ Trim shape for comfort and contoured fit
- √ Can be cold cycle (EtO) sterilized
- Assortment pack (7000D) includes 4 adult, 2 pediatric, 2 infant, and 2 neonate

				All the same
Sensor Type	1 meter	2 meter	3 meter	Extension Cable*
7000A	✓			✓
7000P	✓			✓
7000I	√			✓
7000N	✓			✓



lex System



Neonate Flex System



PureLight

A Rugged Performer



Nonin's Reusable Sensors are perfect for a number of applications where spot checks or short-term monitoring is needed. These durable, quality sensors are designed with precision performance and comfort in mind.

Perfect For Overnight Monitoring

The Flex Sensor Systems are a great solution when extended monitoring is required, such as an overnight study. The Flex sensors are secured with the disposable FlexiWrap® adhesives, providing a comfortable fit without unwanted tightness.



Performance In An Instant



Flexi-Form® II disposable sensors are ideal for emergency situations. Sensors can be easily secured for transport and monitoring on-the-go. EMS personnel particularly like the fact that these sensors perform even in bright light situations.

Sensor Cable Extensions

UNI-RA-O Convenient 6" right angled aďapter



Sensor **Extension Cable** Available in a variety of lengths (1, 3, 6 or 9m)



Please contact your authorized distributor for more information.

Basic Tips for Sensor Selection

Application Criteria	Evaluation Points
Physical Characteristics	General Age of Patient ✓ Infant, Adolescent, Adult Desired Application Site ✓ Finger, Hand, Toe, Foot, Ear, Forehead Patient's Physical Size ✓ Example: A 12 lb. newborn may require an infant sensor rather than a neonatal size
Monitoring Duration	Spot Check ✓ Routine Vitals ✓ Emergency Evaluation Short-Term ✓ Office Procedure ✓ Patient Transport Extended ✓ Therapeutic Assessment ✓ Sleep Study
Amount of Movement	Low Activity Anesthetized, Unconscious Limited Activity Restraints, Tremors General Activity Ambulatory High Activity Exertion Testing
Avoid Cross- Contamination	Disinfect √ Wipe Sensor with Isopropyl Alcohol Sterilize √ Cold Cycle, Ethylene Oxide (EtO)



Your Partner of Choice



With over 18 years of experience and dedication to the design and support of noninvasive monitoring devices, Nonin has helped many medical professionals meet clinical and economic objectives. Trusted by clinicians worldwide, countless tests are performed every day

on thousands of Nonin pulse oximeters in more than 125 countries. Nonin's dedication to technological leadership, precision manufacturing and uncompromised customer support ensures quality products and service you expect.



Nonin Medical, Inc. 13700 1st Avenue North Plymouth, MN 55441-5443, USA Phone: +1.763.553.9968 Toll Free: 800.356.8874 Fax: +1.763.553.7807 info@nonin.com