

PHILIPS

DreamWisp


Nasal mask

Sizing guide

Instructions:

- 1 Print out this Adobe PDF File at 100% on your printer. Use the print accuracy check at right to ensure printing at the proper size.
- 2 Fold the bottom part of this page on the dashed line.
- 3 Cut along the dotted line around the fitting gauge.
- 4 Place the fitting gauge on your nose as shown in the picture below.
- 5 Determine the appropriate size where the outer edge of your nostrils and tip of your nose contact the fitting gauge. (P - Petite, S - Small, M - Medium, L - Large or XL - Extra Large)

Print accuracy check




Tip: For the best performance, use the smallest cushion that fits your nose.

Philips always recommends consulting a physician or registered Sleep Therapist in fitting a mask to ensure a proper sizing and fit.

Fitting gauge

CUT along dotted line



FOLD along dashed line