# StimuScent

## Instructions for use

### Name

StimuScent Smell Training

### Product description

- StimuScent smell training is suitable for training the sense of smell and the brain to (correctly) recognise smells again. Studies show that smell training supports and accelerates recovery of the sense of smell after impairment.
- Each StimuScent smell training set consists of a detailed training schedule and four different, internationally recognised scents.
- People often benefit from smell training, although the end result cannot be predicted in advance.
- Smell training is suitable for adults and for children over 12 years old.
- Insufficient information is available on the use of smell training by pregnant or breastfeeding women. Use during pregnancy or breastfeeding is not recommended.

### Contents

- **Smell training set 1**  
Scent oils (3 ml per bottle).  
Rose, lemon, eucalyptus, clove.
- **Smell training set 2**  
Scent oils (3 ml per bottle).  
Jasmine, mandarin, mint, thyme.
- **Smell training set 3**  
Scent oils (3 ml per bottle).  
Lavender, bergamot, tea tree, rosemary.

### Properties

StimuScent smell training sets each contain 4 different, internationally recognised scents and a detailed training schedule. Each StimuScent Smell training set is a complete set, that makes smell training easy. Smell training can be seen as therapy for the nose. The sense of smell is stimulated and the brain is trained to once again recognise smells correctly.

### Indications/areas of application

Impaired sense of smell.

### How to use (see illustrations)

- 1 Open one of the scent bottles. The order in which the scents are used does not matter.
- 2 Hold the bottle approximately 3 cm under your nose. Do not touch the bottle with your nose.
- 3 Relax. Take one calm, deep breath, inhaling through your nose. This allows the scent to enter your nose effectively. Concentrate on smelling the scent. At the same time, read

the bottle label to see which scent you are smelling. Focus on this particular scent. Do not smell blindly.

- 4 Repeat step 3 four times. You will then have inhaled this smell five times in total.
- 5 Close the bottle. Rest your nose for 5 minutes and move on to the next smell.
- 6 Repeat this process with all the scents in the set. On average, the session will last 30 minutes.
- 7 Repeat the training twice a day (morning/evening), per set, for 12 weeks.

### Frequency and duration of use

- Smell training is preferably performed twice a day, e.g. morning and evening.
- To achieve the best possible results, smell training should be maintained for at least 24 weeks.
- After 12 weeks, switch to a different set of scents. A total of three StimuScent Smell Training sets are available.

### Special attention and precautions

- **Before use, read the hazard statements and precautions on the packaging.**
- **Keep out of reach of children and pets.**
- **Never ingest the scent oil. Swallowing the scent oil causes physical damage.**
- If the scent oil gets on your skin, wash it off with soap and water as soon as possible.
- Do not inhale too quickly, otherwise the scent will not enter your nose properly.
- Concentrate on smelling and take your time. Make sure you know what scent you are smelling and try to be aware of it while smelling. Do not smell blindly.
- Close the bottles tightly after use, otherwise the scent will evaporate.
- Never expose the Smell Training set directly to a heat source.
- Dispose of the bottle and remaining oil at your local municipality's environmental department/waste station.
- Do not use if you are known to be hypersensitive to any of the ingredients.

### Side effects

No side effects have been observed to date.

### Storage

- The expiry date is indicated on the packaging.
- Store at room temperature.
- Can be kept for 3 months after opening.

### Date of last revision of the instructions for use.

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