Air Mate

by **FLAEM**
a healthy lifestyle

GB INSTRUCTIONS FOR USE MANUAL

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E MANUAL DE INSTRUCCIONES

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Mod. P0504EM F400

Aerosol Therapy Device

Congratulations on your purchase and thank you for your confidence in us. Our goal is to fully satisfy our customers by offering them cutting-edge products for treating ailments of the respiratory system. **Please read these instructions carefully and keep them handy for future reference. Use this device only as described in this manual. This is a home medical device for the nebulization and administration of medication prescribed or suggested by your doctor.** Please visit www.flaem.it to see our entire range of Flaem products.

THIS DEVICE COMES EQUIPPED WITH:

- | | |
|---|--|
| <p>Ⓐ - Aerosoltherapy unit (main unit)</p> <ul style="list-style-type: none"> A1 - On/off switch A2 - Air outlet port A3 - Air filter A4 - Nebuliser holder A5 - Carrying handle A6 - Power cord | <p>Ⓒ - Accessories</p> <ul style="list-style-type: none"> C1 - RF8 nebuliser <ul style="list-style-type: none"> C1.1 - Lower part C1.2 - Upper part C2 - Mouthpiece C3 - Adult mask |
| <p>Ⓑ - Connection tube (main unit / nebuliser)</p> | |

⚠ IMPORTANT SAFEGUARDS

- Before using the device for the first time and periodically during the product's life, check the integrity of its structure and of the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorised FLAEM technical assistance centre or to your local dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Keep the cord away from heated surfaces.
- Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.

- Only use the device in dust-free environments; otherwise, therapy may be compromised.
- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to high temperatures.
- Do not put the device near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing in the device.
- Do not obstruct the air vents placed on both sides of the unit.
- During use, always place the unit on a firm surface clear of obstacles.
- Check that there is no material obstructing the air vents before each use.
- Do not place any objects inside the air vents.
- Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.
WARNING: Do not modify this equipment without authorization of the manufacturer.
- The average duration provided for compressor ranges is: F400: 400 hours, F700: 700 hours, F1000: 1000 hours, F2000: 2000 hours.
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

INSTRUCTIONS FOR USE

Before each use, the nebuliser and accessories should be cleaned according to the instructions in the section "CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION". A personal use of the nebuliser and its accessories is recommended in order to avoid any risks of infection.

1. Plug the power cord (A6) into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.
2. Open the nebuliser by turning the upper part (C1.2) counterclockwise.
3. Pour the medication prescribed by the physician into the lower part (C1.1). Close the nebuliser by turning the upper part (C1.2) clockwise.
4. Connect the accessories as shown in the "Assembly diagram".
5. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

6. Start the unit by means of switch (A1) and breathe in and out deeply.

7. When the treatment is finished, turn off and unplug the unit.

WARNING: if after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the compressor air flow; this operation will prevent mold from forming inside the tube.

CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION

Turn off the unit before any cleaning and disconnect the power cord from the mains.

UNIT AND OUTER SURFACE OF THE TUBE

CLEANING

Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).

NEBULIZER AND ACCESSORIES

We recommend a personal use of both the nebulizer and accessories to avoid any risk of infection or contagion. Open the nebulizer by turning the top (C1.2) counterclockwise and disassemble it as shown in the "Connection Diagram" section C1.

SANITIZATION

Before and after each use, sanitize the nebulizer and accessories according to one of the methods described below.

(Method A): The accessories that can be sanitized with the method A are (C1.1-C1.2-C2-C3)

Sanitize the accessories in warm (approximately 40° C) drinking water with mild dishwashing liquid (non abrasive) or in the dishwasher with hot cycle.

(Method B): The accessories that can be sanitized with the method B are (C1.1-C1.2-C2-C3)

Sanitize the accessories by immersion in a solution consisting of 50% water and 50% white vinegar. Finally rinse thoroughly with warm (approximately 40° C) drinking water.

(Method C): The accessories that can be sanitized with the method C are (C1.1-C1.2-C2-C3)

Sanitize the accessories by boiling in water for 20 minutes; it is preferably use use of demineralized or distilled water to avoid scale deposits.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

DISINFECTION

The accessories that can be disinfected are (C1.1-C1.2-C2-C3)

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and

- associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
 - After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
 - Dispose of the solution as recommended by the disinfectant manufacturer.

STERILIZATION

The accessories that can be sterilized are (C1.1-C1.2-C2)

The sterilization procedure described in this paragraph is effective only if all the prescriptions herein are met, all the components involved are previously sanitized and the procedure itself is validated according to ISO 17665-1 standard.

Unit: Steam sterilizer with fractionated vacuum and overpressure according to EN 13060 standard.

Procedure: Package each component to be treated in sterile barrier system or packaging complying with EN 11607 standard. Place the packed components in the steam sterilizer.


Run the sterilization cycle according to the instructions manual of the device, selecting a temperature of 134° C and a duration of 10 minutes.

Storage: Store the sterilized components as per the instructions of the sterile barrier package or system chosen.

AIR FILTERING

The unit is equipped with an air filter (A3) which should be replaced when dirty or when its colour changes. To replace the filter take it out by pulling it as shown in figure. The filter is designed to remain fixed in its housing.

Use original FLAEM accessories only.

SYMBOLS

	Class II device
	Important: check the operating instructions
	“OFF” for part of equipment
	“ON” for part of equipment
	Complies with: European standard EN 10993-1 “Biological Evaluation of Medical Devices” and European Directive 93/42/EEC “Medical Devices”. Phthalates free. In conformity with Reg. (EC) no. 1907/2006
SN	Serial number of device
	Manufacturer
	Type BF applied part
	Alternating current
	Hazard: electrocution. Consequence: Death. Do not use device while taking a bath or shower
	CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments
	TÜV Certification ref. EN 60601-1 3 ^o Ed. - EN60601-1-11
IP21	Enclosure protection rating: IP21. (Protected against solid bodies over 12 mm. Protected against access with a finger; protected against vertically falling water drops.)

AVAILABLE SPARE PARTS

Description	Code
- RF8 Nebuliser +Mouthpiece	ACO597P
- Mask	ACO512P
- Connection tube (1 m)	ACO35
- Air filter (2 pcs)	ACO439P

TECHNICAL FEATURES

Mod. P0504EM F400

Voltage:	230V~ 50Hz 130VA	115V ~ 60Hz	220V ~ 60Hz
Safety certifications:			
Max pressure:	1.8 ± 0.3 bar		
Compressor air output:	9 l/min approx		
Sound level (at 1 m):	54 dB (A) approx		
Operation:	continuous		
Operating conditions:	Temperature: min 10°C; max 40°C RH Air humidity: min 10%; max 95% Atmospheric pressure: min 69KPa; max 106KPa		
Storage conditions:	Temperature: min -25°C; max 70°C RH Air humidity: min 10%; max 95% Atmospheric pressure: min 69KPa; max 106KPa		
Dimensions (W)x(D)x(H):	16x20x12 cm		
Weight:	1.5 Kg		

APPLIED PARTS

Type BF applied parts are: patient accessories (C2,C3)

RF8 NEBULIZER

Medication minimum capacity ml	Medication maximum capacity ml	Operating pressure bar	Delivery ml/min (approx)	MMAD (µm)	BREATHABLE FRACTIONS < 5 µm (FPF) %
2	8	0,8 ⁽¹⁾	0,40 ⁽¹⁾	1,6 ⁽¹⁾⁽²⁾	99,4 ⁽¹⁾⁽²⁾
				4,48 ⁽¹⁾⁽³⁾	55,5 ⁽¹⁾⁽³⁾

⁽¹⁾ Average values measured with 0.9% saline solution in compliance with Flaem internal procedure I29-P07.5.

⁽²⁾ Measured with TSI Aerosizer 3220 computerised laser system in compliance with Flaem internal procedure I29-P07.5.

⁽³⁾ In vitro characterization certified by TÜV Rheinland LGA Product GmbH - Germany in compliance with the new European Standard for aerosol therapy units, Standard EN 13544-1, ANNEX CC. Further details are available on request.

DEVICE DISPOSAL

In conformity with Directive 2012/19/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore deliver (or have delivered) the above waste to an appropriate waste collection centre provided by the local authorities, or deliver it to the dealer when purchasing a new appliance of the same type. Presorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects to the environment and public health due to incorrect waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60601-1-2:2007). Electrical medical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.flaemnuova.it. Flaem reserves the right to make technical and functional changes to the product without notice.

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© 2014 FLAEM NUOVA®
FLAEM NUOVA S.p.A.