



## SmartPeg Kompatibilitätsliste

Implantat Hersteller	System	Durchmesser	Plattform	SmartPeg	Anmerkungen
BIOCONCEPT	BV Tapered	3.5mm	Mini	05	100376
BIOCONCEPT	BV Tapered	4.0mm - 5.0mm	Regular	47	100478
COWELLMEDI	INNO Narrow	3.1mm - 3.3mm	Narrow	28	100432
COWELLMEDI	INNO Internal	3.5mm - 6.0mm	Regular	07	100380
DENTIUM	Simple Line	3.8mm - 4.8mm	Regular	07	100380
MEGAGEN	Any-One	4.5mm - 6.0mm	Regular	07	100380
NEOBIOTECH	IS-Activ	3.5mm-5.0mm	Narrow	05	100376
OSSTEM	TS Serie	3.0mm	Mini	11	100372
OSSTEM	TS Serie	3.5mm	Mini	05	100376
OSSTEM	TS Serie	4.0mm - 5.0mm	Regular	61	100544
OSSTEM	TS Serie	6.0mm - 7.0mm	Ultra Wide	04	100350
HIOSSSEN	ET-Serie	3.5mm	Mini	05	100376
HIOSSSEN	ET-Serie	4.0mm - 6.0mm	Regular	47	100478