



Longest-selling travel CPAP and a proven favorite

Durable, high quality travel CPAP

More power accessories than any other CPAP so you can sleep anywhere

The industry's best warranty (3 year warranty)



PN 400227 REV A

### Transcend 3 miniCPAP Kit

- Transcend 3 miniCPAP Auto
- Padded Transcend travel case with compartments
- Standard six-foot hose
- Multi-plug universal AC power supply with international plug kit
- Mini USB cord
- Transcend 3 miniCPAP Quick Guide



## TRANSCEND<sup>3</sup> minicpap™



	Transcend 3 miniCPAP Auto
Therapy Mode	Auto/Fixed Pressure
EZEX Pressure Relief	YES
Dimensions (cm)	19 x 9.5 x 9.4
Weight (g)	494.4
Sound Pressure Level (dBA)	30
Operating Pressure (cmH <sub>2</sub> O)	4-20
Ramp (0-45 min) Adjustable	YES
Auto Altitude Adjustable (up to 8000 ft)	YES
Adjustable Power Supply Input 100-240VAC, 50-60HZ Power Supply Output 19VDC, 2.6 Amp	YES
AHI and Leak Detect	YES
Compliance Data	YES
TranSync Cloud Compliance	YES

© 2019 Somnetics International, Inc. Protected by U.S. and foreign patents issued and pending. Transcend® 3 miniCPAP™ is a registered trademark of Somnetics International, Inc.

*new*  
TRANSCEND<sup>3</sup> minicpap

Freedom to sleep anywhere™



[www.MyTranscend.com](http://www.MyTranscend.com)

# THE miniCPAP™ JUST GOT BETTER.



The original Transcend miniCPAP™ has been a top-selling travel CPAP for years and now it's even better with the latest upgrades. The reimagined Transcend 3 miniCPAP™ maintains the small, lightweight build you've come to expect in a miniCPAP with upgraded features like a modern design, a flat silicone base for better stability, and a swivel nozzle for a better CPAP experience.

## Mini Size & Weight

About the size of a soda can, and weighing just a pound, the Transcend 3 miniCPAP is easy to pack and carry on your next trip. It's even FAA-approved so you can sleep comfortably on your next flight. And you can add a battery pack that's about the size of a deck of cards so you can sleep anywhere.



## Durable Quality Design

The Transcend 3 miniCPAP is designed for travel. This robust miniCPAP supports your active lifestyle so you have the freedom to travel AND sleep well no matter where you go. It is even backed by an industry-leading warranty.

## The Most Power Options

Transcend 3 miniCPAP has power options that make it easy for you to sleep anywhere. From five star hotels to the rugged outdoors, we've got you covered. Choose from wide variety of power accessories that offer unmatched power options to support your portable CPAP therapy on the go.

## GET THE MOST OUT OF YOUR TRANSCEND 3 MINICPAP WITH ACCESSORIES

### P8 Battery

Small and lightweight, it's perfect for sleeping with your Transcend 3 miniCPAP outdoors or on your next flight. This FAA-approved accessory ensures you always have the power you need.

MOST  
POPULAR

### DC Mobile Power Adaptor

Recharge the P8 Battery or power your Transcend 3 miniCPAP from any 12-volt outlet.



### Solar Battery Charger

Perfect for camping and remote getaways. This foldable lightweight solar panel recharges your P8 battery with the sun so you can sleep under the stars.



TRANSCEND<sup>3</sup>  
minicpap™

