



# Contec Filtered Alcohols

70% IPA or Denatured Ethanol blended with purified water

Contec alcohols are a blend of 70% v/v Isopropanol or denatured ethanol with 30% purified water. The alcohol blend is 0.2 micron filtered, filled and bagged in a Grade C (ISO Class 7) cleanroom.

Contec’s ethanol is only denatured with IPA creating a very low residue product.

Supplied as either trigger sprays fitted with a protected system, which ensures sterility throughout use, or 5L capped containers for when larger volumes are required.



FEATURES	BENEFITS
GMP manufactured under Grade A air flow in a Grade C cleanroom	<ul style="list-style-type: none"> <li>Ensures the alcohol, container and packaging are free from contamination and particulates</li> </ul>
No denaturants other than IPA	<ul style="list-style-type: none"> <li>Completely residue free</li> </ul>
Trigger sprayer and “bag in bottle” protected system	<ul style="list-style-type: none"> <li>Prevents bottle contents from becoming contaminated during use</li> <li>No need to discard unused product during shelf life</li> <li>Bottles can be completely emptied eliminating wastage</li> </ul>
Trigger sprayer can be set to jet or spray	<ul style="list-style-type: none"> <li>Large droplet size reduces the risk of inhalation and provides good surface coverage</li> </ul>
Double bagged packed in linear tear packaging	<ul style="list-style-type: none"> <li>Each bag is easy to open even when wearing gloves</li> <li>Facilitates transfer disinfection into cleanroom</li> </ul>

## Product Specification

70% v/v IPA or denatured ethanol with 30% purified water

Ready to use with a 3 year shelf life

0.2 micron filtered

Double bagged in polyethylene linear tear bags

Tamper evident caps and triggers

Certificates of analysis provided with every batch.

Full supporting documentation including efficacy information (EN1276 and EN1650) is available on request.



Copyright © 2016 Contec, Inc. All rights reserved.



**PROCLEANROOM**  
BUY YOUR SUPPLIES ONLINE

## Product Application

Ready to use so can be applied directly to the surface ensuring even coverage. Physically remove contaminants from the surface using an appropriate cleanroom wipe and a recommended wiping technique for optimum contamination control.

Part No.	Description	Size	Packaging
FBT170I	Non-sterile 70% isopropanol in purified water with trigger sprayer	1L	6/case
FBC570I	Non-sterile 70% isopropanol in purified water	5L	2/case
FBT170DE	Non-sterile 70% denatured ethanol in purified water with trigger sprayer	1L	6/case
FBC570DE	Non-sterile 70% denatured ethanol in purified water	5L	2/case