

STERILE NEUTRAL DETERGENT WIPES

**IN WFI 68GSM PRE-IMPREGNATED – READY TO USE
ITEM 192900 (10 X 10 WIPES)**

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANT APPLICATIONS.

AGMA HAVE CONSIDERABLE EXPERIENCE IN SURFACTANT TECHNOLOGY AND STERILE CLEANING / DISINFECTION. WE HAVE COMBINED THIS EXPERTISE TO DEVELOP A STERILE NEUTRAL DETERGENT SPECIFICALLY FOR USE IN CLEANROOMS AND PHARMACIES.

WITH HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES, AGMA STERILE NEUTRAL DETERGENT IS IDEAL WHEN A DEEP CLEAN IS REQUIRED. FORMULATED TO BE LOW FOAMING IT DOES NOT LEAVE SIGNIFICANT RESIDUES OR SMEARS. IN ADDITION, ITS PH HELPS PROTECT MORE DELICATE WORK SURFACES AND THOSE PRONE TO CORROSION.

- **MADE WITH WFI QUALITY WATER.**
- **0.2 MICRON FILTERED.**
- **GAMMA IRRADIATED FOR STERILITY.**
- **LOW FOAMING.**
- **HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES.**
- **MANUFACTURED IN CLASS C/D CLEANROOMS.**
- **REMOVES GREASES, DISINFECTANT RESIDUE ETC.**
- **DOES NOT LEAVE SIGNIFICANT RESIDUES.**
- **68GSM, LOW LINTING POLYCELLULOSE WIPE.**
- **28CM X 38CM FOR LARGE SURFACE AREA COVERAGE.**
- **FULLY TRACEABLE WITH SUPPORTING ANALYSIS.**
- **2 YEAR SHELF LIFE.**



DIRECTION FOR USE

AGMA NEUTRAL DETERGENT WIPES ARE MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1. IDEAL FOR USE ON CLEANROOM HARD SURFACES TO REMOVE RESIDUES AND FOR USE BETWEEN ROTATING OF DISINFECTANTS. AN INITIAL DEEP CLEAN SHOULD BE PERFORMED PRIOR TO THE USE OF ANY DISINFECTANT TO HELP THE REMOVAL OF ANY SOILS OR HEAVY RESIDUE BUILD UP.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-02/25