

IPA 70% IN WFI 4YPINS

ITEM 1213 (10 X 50 WIPES)

0.2 MICRON FILTERED

IPA 70% IN WFI IS A PROVEN MEANS FOR WIPING INTO PHARMACEUTICAL & HOSPITAL CLEANROOMS.

READY TO USE WIPES ARE A SIMPLE MEANS OF DISINFECTING A VARIETY OF SURFACES AND CLEANROOM ITEMS. THEY HAVE THE ADDED BENEFIT OF BEING ABLE TO PHYSICALLY REMOVE PARTICLES AND SPORES IN ADDITION TO PROVIDING GOOD DISINFECTION AGAINST VEGETATIVE MICROBES.

COST IN USE, TAKING INTO ACCOUNT OPERATOR TIME, GIVES FINANCIAL BENEFITS WHEN COMPARED TO DRY WIPES AND SPRAYS.

- READY TO USE-SAVES TIME.
- STANDARD QUALITY PACKS OF 50 WIPES.
- 0.2 MICRON FILTERED.
- RE-SEALABLE TO PREVENT DRYING OUT.
- COST EFFECTIVE IN USE.
- BROAD SPECTRUM ACTIVITY.
- MADE WITH WFI QUALITY WATER.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- POLYPROPYLENE Z FOLDED, 30GSM
25CM X 28CM WIPES.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- 2 YEAR SHELF LIFE.



DIRECTION FOR USE:

AGMA DE 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS

ISSUE NO. 1-03/25

STERILE IPA 70% IN WFI - 4YPINS

30GSM PRE-IMPREGNATED – READY TO USE / EASY TO USE.

ITEM 31824 (10 X 50 WIPES)

STERILE IPA 70% IN WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

READY TO USE WIPES ARE A SIMPLE MEANS OF DISINFECTING A VARIETY OF SURFACES AND CLEANROOM ITEMS. THEY HAVE THE ADDED BENEFIT OF BEING ABLE TO PHYSICALLY REMOVE PARTICLES AND SPORES IN ADDITION TO PROVIDING GOOD DISINFECTION AGAINST VEGETATIVE MICROBES.

COST IN USE, TAKING INTO ACCOUNT OPERATOR TIME, GIVES FINANCIAL BENEFITS WHEN COMPARED TO DRY WIPES AND SPRAYS.

- READY TO USE - SAVES TIME.
- POLYPROPYLENE BASED WIPE.
- 0.2 MICRON FILTERED.
- 25CM X 25CM - 30GSM WIPES.
- 10 X 50 WIPES.
- RE-SEALABLE TO PREVENT DRYING OUT.
- WFI QUALITY WATER – HIGH PURITY WATER.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATED FOR STERILITY.
- 2 YEAR SHELF LIFE.



DIRECTION FOR USE:

AGMA IPA 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS, AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE THE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-03/25