



# WiZARD 310/320




Nasal Mask / Full Face Mask  
Sizing Gauge



WiZARD 310  
Nasal Mask



WiZARD 320  
Full Face Mask



1. Please print out this sizing gauge and ensure it is printed in scale.
2. Find your mask model and cut along dotted red line.
3. Place the fitting gauge on your nose refer to the left picture and determine appropriate cushion size.

**APEX**

0(inch)

1

0(cm)

1

2

This ruler is for scaling purposes.

Please print out this sizing gauge at 100%.