

















**OXYCLOG-EBL** 



**OXYCLOG-BLU** 



**OXYCLOG-EGN** 



**OXYCLOG-FUX** 







**OXYCLOG-LGN** 



**OXYCLOG-LIC** 



**OXYCLOG-RED** 



**OXYCLOG-WHT** 

The model OXYCLOG is a clog for the operation room with an adjustable strap to guarantee a perfect fit. The anatomic shaped foot bed gives optimal support.

This shoe easily meets every requirement of the OR. This clog can be chemically or UV sterilized and is autoclavable. The shoe is also wachine washable at 90°C or 194°F.

The slip resistant outsole provides a good grip both on wet and dry surfaces. The closed top is to protect the foot, but lateral vents take care of the needed ventillation.

Oxypas is part of the Cortina Group

Cortina NV - Meersbloem Melden 42 - 9700 Oudenaarde, Belgium

E-mail: contact@oxypas.com

Phone: (32) 55 33 44 44 - Fax: (32) 55 33 56 46

Website: www.oxypas.com

# SOXYPAS Oxyclog

#### CE standard

CE EN 20347 certified.



### Autoclavable

Till 135° C or 275°F.



# Washable at 90°C or



# Chemical and UV sterilizable



#### Shock absorbing

The construction of the heel makes it very shock absorbing.



#### Thermoplastic elastomer

The shoe is made of TPE. This is why the shoe remains flexible and can be sterilised on high temperatures without any deformation.





#### Weight: 270 gr

The average weight of the shoe is 9.52 ounces or 270 grams.



#### Anti-static ESD

The shoe has a plug that provides dissipation of static electricity to be confrom with the ESD standard.



#### Slip resistant SRA

The outsole provides excellent grip, both on wet and dry floors and complies with the SRA standard.

## Recommended for

- Operating room personnel
- Professionals who help and assist patients in the shower
- Healthcare professionals
- Nurses and caregivers
- Cleaning personnel
- Laundry personnel

## Sizes

EU	35/36	37/38	39/40	41/42	43/44*	45/46*
UK	2.5/3.5	4/5	5.5/6.5	7/8	9/9.5	10.5/11
US	3/4	4.5/5.5	6.5/7.5	8/9	10/10.5	11.5/12

(\*) The colors fuchsia and purple are only available from size 36 till 42